
REMOVABLE "Y" LUER LOCK ADAPTOR INSTRUCTIONS FOR USE

INTENDED USE: To provide dual access to a single lumen hemodialysis catheter.

INSTALLATION PROCEDURE:

1. Install "Y" Adaptor male luer onto female luer of catheter.
2. Turn base of adaptor clockwise until secure.
3. Flush adaptor with saline to ensure connection is secure.
4. Connect to bloodlines to begin treatment.

POTENTIAL COMPLICATIONS:

- Blood Loss
- Infections

WARNINGS & PRECAUTIONS:


- To prevent accidents, assure the security of all caps and bloodline connections prior to and between treatments.
- Discard biohazard according to facility protocol.
- Examine catheter female luer after each treatment for damage.
- Federal (USA) law restricts this device to sale by or on the order of a physician.
- Product is sterile in unopened, undamaged package. STERILIZED BY ETHYLENE OXIDE

STERILE	EO
---------	----
- Single use only. DO NOT RE-USE. Re-Use may lead to infection or illness/injury. 
- Do not use if package is damaged. 
- Do not resterilize. 
- Caution when using this device. Be aware of sharps.

- Do not overtighten. Do not proceed if resistance is felt or interaction between components is failing.
- Examine the device and connections before each treatment.

Medcomp® does not recommend a particular technique for the use of this device. The physician should evaluate the appropriateness of the device according to the individual patient conditions and his or her medical training and experience.

SYMBOL TABLE

5.1.1		Manufacturer *
5.3.4		Keep Dry *
5.6.3		Non-pyrogenic *
5.3.2		Keep Away from Sunlight *
5.2.3		Sterilized Using Ethylene Oxide *
5.2.8		Do Not Use if Package is Damaged *
5.1.4		Use-by Date *
5.2.6		Do Not Resterilize *
5.1.5		Batch/Lot Number *
5.1.6		Catalogue Number *
5.4.4		Caution, consult Accompanying Documents *
5.4.2		Do Not Re-use *
	Rx Only	Prescription Use Only ***

***This symbol is in accordance with ISO 15223-1.**

***** FDA guidance Use of Symbols in Labeling.**

WARRANTY

Medcomp® WARRANTS THAT THIS PRODUCT WAS MANUFACTURED ACCORDING TO APPLICABLE STANDARDS AND SPECIFICATIONS. PATIENT CONDITION, CLINICAL TREATMENT, AND PRODUCT MAINTENANCE MAY AFFECT THE PERFORMANCE OF THIS PRODUCT. USE OF THIS PRODUCT SHOULD BE IN ACCORDANCE WITH THE INSTRUCTIONS PROVIDED AND AS DIRECTED BY THE PRESCRIBING PHYSICIAN.

Because of continuing product improvement, prices, specifications, and model availability are subject to change without notice. Medcomp® reserves the right to modify its products or contents without notice.

Medcomp® is a trademark of Medical Components, Inc. registered in the United States

 **Medical Components, Inc.**
 1499 Delp Drive
 Harleysville, PA 19438 U.S.A.
 Tel: 215-256-4201
 Fax: 215-256-1787
 www.medcompnet.com

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Medcomp® est une marque de Medical Components, Inc. enregistrée aux États-Unis.