medCOMP®

END CAP

INSTRUCTIONS FOR USE

INTENDED USE:

- To keep clean and protect catheter luer between treatments.

- Allowing access to the catheter.

DESCRIPTION:

- Self-sealing diaphragm under a plastic cap with a luer connector.

POTENTIAL COMPLICATIONS:

- Accidental disconnection
- Bleeding
- Air embolism
- Infection

WARNINGS:

- Federal Law (USA) restricts the device to sale by or on the order of a physician.

- This device is for single use only. 

- Do not resterilize by any method. 

- Re-use may lead to infection or illness/injury.

- The manufacturer shall not be liable for any damages caused by re-use or resterilization of the device.

- Contents sterile and non-pyrogenic in un-opened , undamaged package.

- STERILIZED BY ETHYLENE OXIDE

STERILE	EO
---------	----

- Do not use if package is opened or damaged. 

- Do not use if any sign of product damage is visible.

INSTRUCTIONS FOR USE:

- After locking the catheter with locking solution.

- Remove locking solution syringe.

- Remove dust cap, if applicable.

- Thread end cap onto catheter luer.

- Use a slight twisting motion to ensure that the end cap is fully seated/connected onto the luer.

PRECAUTIONS:

- Do not overtighten - may cause cracking of the luer and make removal difficult.

- Do not use hemostats or another instrument to remove cap - may cause cracking.

- Use only luer lock connectors.

- Discard biohazard according to facility protocol.

Medcomp® does not recommend a particular technique for the use of this device. The physician should evaluate their appropriateness according to individual patient conditions and his or her medical training and experience.

WARRANTY

Medcomp® WARRANTS THAT THIS PRODUCT WAS MANUFACTURED ACCORDING TO APPLICABLE STANDARDS AND SPECIFICATIONS. PATIENT CONDITION, CLINICAL TREATMENT, AND PRODUCT MAINTENANCE MAY EFFECT THE PERFORMANCE OF THIS PRODUCT. USE OF THIS PRODUCT SHOULD BE IN ACCORDANCE WITH THE INSTRUCTIONS PROVIDED AND AS DIRECTED BY THE PRESCRIBING PHYSICIAN.

Because of continuing product improvement, prices, specifications, and model availability are subject to change without notice. Medcomp® reserves the right to modify its products or contents without notice.

Medcomp® is a registered trademark of Medical Components, Inc.

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Medcomp® jest registrirani zaštitni znak društva Medical Components, Inc.

SYMBOL TABLE

5.1.1 	Manufacturer *
5.3.4 	Keep Dry *
5.4.2 	Do Not Re-use *
5.6.3 	Non-pyrogenic *
5.3.2 	Keep Away from Sunlight *
5.2.3 	Sterilized Using Ethylene Oxide *
5.2.8 	Do Not Use if Package is Damaged *
5.1.4 	Use-by Date *
5.2.6 	Do Not Resterilize *
5.1.5 	Batch/Lot Number *
5.1.6 REF	Catalogue Number *
5.1.2 	Authorized Representative in the European Community *
Rx Only	Prescription Use Only ***

*** This symbol is in accordance with ISO 15223-1.**

***** FDA guidance Use of Symbols in Labeling.**

EC REP

MPS Medical Product Service GmbH
 Borngasse 20
 35619 Braunfels
 Germany

Medical Components, Inc.
 1499 Delp Drive
 Harleysville, PA 19438 U.S.A.
 Tel:215-256-4201
 Fax:215-256-1787
 www.medcompnet.com