4.0% DuraLock-C™  
CATHETER LOCK SOLUTION

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## **4.0% DuraLock-C™**

### **PRODUCT DESCRIPTION:**

DuraLock-C™ is supplied as a clear, plastic, prefilled syringe, packaged in a pouch. Each DuraLock-C™ prefilled syringe contains a sterile, clear and colorless solution containing sodium citrate; pH is adjusted with citric acid.

### **INDICATIONS AND CLINICAL USE:**

DuraLock-C™ is indicated for use in maintaining patency of Hemodialysis Catheters.

### **ACTION:**

Trisodium Citrate causes anticoagulation by chelating ionized Calcium into a soluble complex. Calcium is an integral ion involved in the clotting cascade. Local removal by Citrate prevents the activation of clotting cofactors, factor X and prothrombin and the ultimate formation of fibrin. Systemic anticoagulation does not occur.

### **INTENDED USE:**

Short-Term, between 24hours-72hours, catheter lock solution.

### **WARNING:**

If catheter is not sufficiently locked, clots can form and thrombosis may occur.

### **CONTRAINDICATIONS:**

Patients with known allergies or hypersensitivity to constituents.

### **STORAGE:**

Duralock-C™ should be stored between 15-30°C and away from exposure to direct sunlight and moisture.

## PRECAUTIONS:

- Injection volume in excess of catheter manufacture's specific lumen priming volume should not be used.
- For use with adult patients only.
- Use aseptic technique.
- Do not apply by direct intravenous injection or added to an infusion.
- Only to be used when exact catheter lumen volumes are known.
- Single use only, do not re-use. Re-use of single use devices creates a potential risk of contamination. Ⓢ This may lead to patient infection, which may cause illness or death of the patient or user.
- For prescription use only and to be used by a skilled or trained professional only. **Rx Only**
- Do not use if package is opened or damaged. Ⓢ
- Do not use if the cap on the syringe is not intact or damaged.
- Before using, ensure that catheter lumen, extensions and luer do not contain cracks, swelling or other signs of damage before and after each treatment session.
- Prior to use, visually inspect the contents of each prefilled syringe for clarity, particulate matter, precipitate, discoloration and leakage. **NOTE:** Solutions showing haziness, particulate matter, precipitate, discoloration or leakage should not be used.
- Contents are sterile and non-pyrogenic in unopened, undamaged package. Sterilized by gamma irradiation. **STERILE R** Ⓢ
- Lot number and expiration date are located on pouch and syringe label.
- If skin contact occurs, clean per facility's procedures.

## ADVERSE REACTIONS:

No known adverse reactions when the product is used as indicated. Paraesthesia and/or dysgeusia may occur if the product unintentionally passes into the vein.

## INSTRUCTIONS FOR USE (Use Aseptic Technique):

1. Prior to initiation of a venous access session, **aspirate and discard the solution that was instilled in the Hemodialysis Catheter during the previous session**
2. At the completion of each venous access session, flush the lumens of the catheter with 5 mL of sterile 0.9% Sodium Chloride solution to remove any blood, medication or other substance remaining in the catheter.
3. Aseptically open package containing DuraLock-C™ prefilled syringe(s). Inspect the contents of the syringe(s) according to PRECAUTIONS.
4. Hold the syringe with the cap end up and carefully twist the cap off (Fig. 1).
5. To remove any air bubbles, point luer end of syringe up and gently tap the syringe with fingers until air bubbles rise to top of the syringe.
6. Keep luer end of syringe pointing up and slowly push the plunger up to force air bubbles out of the syringe.
7. **Properly identify priming volume** and with luer end of syringe pointing up, slowly **push plunger to the line on the syringe that matches the priming volume** (Prescribed by the catheter manufacturer) (Fig. 2, Fig.3).
8. Connect the syringe to the catheter. Use only with compatible luer lock connectors (Fig. 4).
9. Slowly, over a 5 to 10 second period, **instill the exact priming volume**, as specified by the catheter manufacturer, of DuraLock-C™ solution into each catheter lumen, as a locking agent.
10. Remove and discard unused portions and empty syringes according to the institution's biohazardous waste policy.


Fig. 1


Fig. 2


Fig. 3


Fig. 4

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# Medical Components, Inc.
*Harleysville, PA 19438 U.S.A.*