




CATHETER CLAMP SLEEVE



1. After tunneling, slide clamp sleeve over catheter lumen.




2. Install clamp over clamp sleeve.




3. Install adapter according to supplied *Instructions for Use*.

FOR BEST RESULTS, MOVE POSITION OF CLAMP PERIODICALLY.


CATHETER CLAMP SLEEVE



1. After tunneling, slide clamp sleeve over catheter lumen.







2. Install clamp over clamp sleeve.



3. Install adapter according to supplied *Instructions for Use*.

FOR BEST RESULTS, MOVE POSITION OF CLAMP PERIODICALLY.

WARNINGS:

- Federal Law (USA) restricts the device to sale by or on the order of a physician.
- Product is sterile in unopened, undamaged package. Sterilized by ethylene oxide. 
- Single use only. DO NOT RE-USE. Re-use may lead to infection or illness/injury. 
- Do not use if package is damaged. 
- Do not resterilize. 














Note: Discard biohazard according to facility protocol.

WARRANTY

Medcomp® WARRANTS THAT THIS PRODUCT WAS MANUFACTURED ACCORDING TO APPLICABLE STANDARDS AND SPECIFICATIONS, PATIENT CONDITION, CLINICAL TREATMENT, AND PRODUCT MAINTENANCE MAY AFFECT THE PERFORMANCE OF THIS PRODUCT. USE OF THIS PRODUCT SHOULD BE IN ACCORDANCE WITH THE INSTRUCTIONS PROVIDED AND AS DIRECTED BY THE PRESCRIBING PHYSICIAN.

Because of continuing product improvement, prices, specifications, and model availability are subject to change without notice. Medcomp® reserves the right to modify its products or contents without notice.





Medcomp® is a registered trademark of Medical Components, Inc.

5.1.1 	Manufacturer *	5.6.3 	Non-pyrogenic *
5.1.4 	Use-by Date *	5.3.2 	Keep Away from Sunlight *
5.4.2 	Do Not Re-use *	5.2.6 	Do Not Resterilize *
5.3.4 	Keep Dry *	Rx Only	Prescription Use Only ***
5.1.5 	Batch/Lot Number *	5.1.6 REF	Catalogue Number *
5.2.8 	Do Not Use if Package is Damaged *		
5.1.2 	Authorized Representative in the European Community *		
5.2.3 	Sterilized Using Ethylene Oxide *		
5.4.3 	Consult Instructions for Use *		
5.4.4 	Caution, consult Accompanying Documents		

*This symbol is in accordance with ISO 15223-1.

*** FDA guidance Use of Symbols in Labeling.

WARNINGS:

- Federal Law (USA) restricts the device to sale by or on the order of a physician.
- Product is sterile in unopened, undamaged package. Sterilized by ethylene oxide. 
- Single use only. DO NOT RE-USE. Re-use may lead to infection or illness/injury. 
- Do not use if package is damaged. 
- Do not resterilize. 














Note: Discard biohazard according to facility protocol.

WARRANTY

Medcomp® WARRANTS THAT THIS PRODUCT WAS MANUFACTURED ACCORDING TO APPLICABLE STANDARDS AND SPECIFICATIONS, PATIENT CONDITION, CLINICAL TREATMENT, AND PRODUCT MAINTENANCE MAY AFFECT THE PERFORMANCE OF THIS PRODUCT. USE OF THIS PRODUCT SHOULD BE IN ACCORDANCE WITH THE INSTRUCTIONS PROVIDED AND AS DIRECTED BY THE PRESCRIBING PHYSICIAN.

Because of continuing product improvement, prices, specifications, and model availability are subject to change without notice. Medcomp® reserves the right to modify its products or contents without notice.

Medcomp® is a registered trademark of Medical Components, Inc.

5.1.1 	Manufacturer *	5.6.3 	Non-pyrogenic *
5.1.4 	Use-by Date *	5.3.2 	Keep Away from Sunlight *
5.4.2 	Do Not Re-use *	5.2.6 	Do Not Resterilize *
5.3.4 	Keep Dry *	Rx Only	Prescription Use Only ***
5.1.5 	Batch/Lot Number *	5.1.6 REF	Catalogue Number *
5.2.8 	Do Not Use if Package is Damaged *		
5.1.2 	Authorized Representative in the European Community *		
5.2.3 	Sterilized Using Ethylene Oxide *		
5.4.3 	Consult Instructions for Use *		
5.4.4 	Caution, consult Accompanying Documents		


*This symbol is in accordance with ISO 15223-1.

*** FDA guidance Use of Symbols in Labeling.

 **Medical Components, Inc.**
 1499 Delp Drive
 Harleysville, PA 19438 U.S.A.
 Tel:215-256-4201
 Fax:215-256-1787
 www.medcompnet.com


 MPS Medical Product Service GmbH
 Borngasse 20
 35619 Braunfels
 Germany




 MPS Medical Product Service GmbH
 Borngasse 20
 35619 Braunfels
 Germany



 **Medical Components, Inc.**
 1499 Delp Drive
 Harleysville, PA 19438 U.S.A.
 Tel:215-256-4201
 Fax:215-256-1787
 www.medcompnet.com