




Silicone Material	Step Tip Tip Design	8F, 12.5F French Sizes	Straight, Pre-Curved Configurations	Antegrade Insertion Type
-----------------------------	-------------------------------	----------------------------------	---------------------------------------------------	------------------------------------



CIRCLE "C" INTERNAL LUMEN DESIGN



LONGER LENGTHS
Enable adequate treatment with tip placement in the IVC via the femoral approach.



SILICONE MATERIAL
Provides biocompatibility, kink resistance and improved patient comfort.



ANGLED ARTERIAL TIP
Eases insertion and reduces vessel trauma.

BASIC SET		
SL18P	8F X 18CM HEMO-CATH® LT, STRAIGHT	5/BOX
SL24P	8F X 24CM HEMO-CATH® LT, STRAIGHT	5/BOX

CATHETER ONLY		
SL28C	12.5F X 28CM HEMO-CATH® LT, STRAIGHT	10/BOX
SL32C	12.5F X 32CM HEMO-CATH® LT, STRAIGHT	10/BOX
SL28CPC	12.5F X 28CM HEMO-CATH® LT, PRE-CURVED	10/BOX
SL32CPC	12.5F X 32CM HEMO-CATH® LT, PRE-CURVED	10/BOX

BASIC SET		
MC101241	12.5F X 15CM HEMO-CATH® LT, STRAIGHT	5/BOX
MC101242	12.5F X 18CM HEMO-CATH® LT, STRAIGHT	5/BOX
MC101243	12.5F X 24CM HEMO-CATH® LT, STRAIGHT	5/BOX
SL28	12.5F X 28CM HEMO-CATH® LT, STRAIGHT	5/BOX
SL32	12.5F X 32CM HEMO-CATH® LT, STRAIGHT	5/BOX
SL28PC	12.5F X 28CM HEMO-CATH® LT, PRE-CURVED	5/BOX
SL32PC	12.5F X 32CM HEMO-CATH® LT, PRE-CURVED	5/BOX

FULL TRAY		
MCDLTSL28	12.5F X 28CM HEMO-CATH® LT, STRAIGHT	5/BOX
MCDLTSL32	12.5F X 32CM HEMO-CATH® LT, STRAIGHT	5/BOX
MCDLTSL28PC	12.5F X 28CM HEMO-CATH® LT, PRE-CURVED	5/BOX
MCDLTSL32PC	12.5F X 32CM HEMO-CATH® LT, PRE-CURVED	5/BOX

8F BASIC SET CONTAINS:

- (1) Catheter
- (1) 10cc Syringe
- (1) Guidewire
- (1) Adhesive Wound Dressing
- (1) Hemo-Cath® Hub Clip
- (1) Introducer Needle
- (2) Injection Caps
- (1) Trocar
- (1) Peelable Sheath Introducer
- (1) Patient I.D. Card

CATHETER ONLY CONTAINS:

- (1) Catheter
- (2) Injection Caps
- (1) Hemo-Cath® Hub Clip
- (1) Patient I.D. Card

BASIC SET CONTAINS:

- (1) Catheter
- (2) Injection Caps
- (1) Scalpel
- (1) Peelable Sheath Introducer
- (1) Introducer Needle
- (1) Guidewire
- (1) Adhesive Wound Dressing
- (1) Hemo-Cath® Hub Clip
- (1) Trocar
- (1) Patient I.D. Card

FULL TRAY CONTAINS:

- (1) Catheter
- (1) Dilator
- (1) Peelable Sheath Introducer
- (4) 4" x 4" Gauze
- (1) Chlorhexidine
- (2) 5cc Syringes
- (3) Safety Hypodermic Needles
- (1) 5cc Lidocaine
- (1) Safety Introducer Needle
- (1) Surgical Gloves
- (2) 10cc Syringes
- (1) Safety Scalpel
- (1) Guidewire
- (1) 2-0 Silk Suture
- (1) Fenestrated Drape
- (1) Hemostat
- (1) Adhesive Wound Dressing
- (2) Injection Caps
- (1) Trocar
- (1) Hemo-Cath® Hub Clip
- (1) Patient I.D. Card