




# LT SILICONE CVC

<b>Single, Double, Triple</b> Lumens	<b>90cm</b> Length	<b>6F, 7F, 9.5F, 11F, 12.5F, 14F</b> French Sizes
---	-----------------------	--



**SOFT, SILICONE MATERIAL**  
Provides biocompatibility and kink resistance, enhancing long term performance.

**POLYESTER CUFF**  
Promotes tissue in-growth, securing the catheter and reducing infection.

	BASIC SET	
MR28010611	6F X 90CM SINGLE LUMEN LT CVC	10/BOX
MR28010711	7F X 90CM SINGLE LUMEN LT CVC	10/BOX
MR28010911	9.5F X 90CM SINGLE LUMEN LT CVC	10/BOX
MR28010721	7F X 90CM DOUBLE LUMEN LT CVC	10/BOX
MR28010921	9.5F X 90CM DOUBLE LUMEN LT CVC	10/BOX
MR28011121	11F X 90CM DOUBLE LUMEN LT CVC	10/BOX
MR28011421	14F X 90CM DOUBLE LUMEN LT CVC	10/BOX
MR28011231	12.5F X 90CM TRIPLE LUMEN LT CVC	10/BOX

**BASIC SET CONTAINS:**  
 (1) Infusion Catheter  
 (1) Peelable Introducer  
 (1) Tunneler  
 (1) Mini Scalpel  
 (1) Guidewire  
 (1|2|3) Injection Cap(s)  
 (1) Introducer Needle