


MULTI-CATH

Single, Double,
Triple, Quad
Lumens

5cm, 8cm, 13cm,
15cm, 20cm, 30cm
Lengths

14GA, 16GA, 18GA,
20GA, 24GA, 4F, 5F,
5.5F, 7F, 8F, 8.5F
Sizes



SOFT, FLEXIBLE, RADIOPAQUE TIPS
Less traumatic to the vessel and allow for fluoroscopic confirmation of catheter placement.

PRINTED DEPTH MARKINGS
Assist in accurate central venous placement.

MULTIPLE CONFIGURATIONS
Available in a variety of lengths, French sizes, and lumen configurations.

	SINGLE LUMEN BASIC SET	
CSP-1406	14GA X 15CM MULTI-CATH, SINGLE LUMEN	20/BOX
CSP-1408	14GA X 20CM MULTI-CATH, SINGLE LUMEN	20/BOX
CSP-1412	14GA X 30CM MULTI-CATH, SINGLE LUMEN	20/BOX
CSP-1606	16GA X 15CM MULTI-CATH, SINGLE LUMEN	20/BOX
CSP-1608	16GA X 20CM MULTI-CATH, SINGLE LUMEN	20/BOX
CSP-1612	16GA X 30CM MULTI-CATH, SINGLE LUMEN	20/BOX
CSP-1806	18GA X 15CM MULTI-CATH, SINGLE LUMEN	20/BOX
CSP-1808	18GA X 20CM MULTI-CATH, SINGLE LUMEN	20/BOX
CSP-1812	18GA X 30CM MULTI-CATH, SINGLE LUMEN	20/BOX
CSP-2005	20GA X 13CM MULTI-CATH, SINGLE LUMEN	20/BOX
CSP-2405	24GA X 13CM MULTI-CATH, SINGLE LUMEN	20/BOX

	DOUBLE LUMEN BASIC SET	
MC2LP-35S	4F X 5CM MULTI-CATH, DOUBLE LUMEN	10/BOX
MC2LP-38S	4F X 8CM MULTI-CATH, DOUBLE LUMEN	10/BOX
MC2LP-313S	4F X 13CM MULTI-CATH, DOUBLE LUMEN	10/BOX
MC2LP-330S	4F X 30CM MULTI-CATH, DOUBLE LUMEN	10/BOX
MC2LP-58S	5F X 8CM MULTI-CATH, DOUBLE LUMEN	10/BOX
MC2LP-5S	5F X 13CM MULTI-CATH, DOUBLE LUMEN	10/BOX
MC2LP-530S	5F X 30CM MULTI-CATH, DOUBLE LUMEN	10/BOX
MC2L-6S	7F X 15CM MULTI-CATH, DOUBLE LUMEN	10/BOX
MC2L-8S	7F X 20CM MULTI-CATH, DOUBLE LUMEN	10/BOX
MC2L-12S	7F X 30CM MULTI-CATH, DOUBLE LUMEN	10/BOX

	TRIPLE LUMEN BASIC SET	
MC3LP-58S	5.5F X 8CM MULTI-CATH, TRIPLE LUMEN	10/BOX
MC3LP-5S	5.5F X 13CM MULTI-CATH, TRIPLE LUMEN	10/BOX
MC3LP-530S	5.5F X 30CM MULTI-CATH, TRIPLE LUMEN	10/BOX
MC3L-6S	7F X 15CM MULTI-CATH, TRIPLE LUMEN	10/BOX
MC3L-8S	7F X 20CM MULTI-CATH, TRIPLE LUMEN	10/BOX
MC3L-12S	7F X 30CM MULTI-CATH, TRIPLE LUMEN	10/BOX

	QUAD LUMEN BASIC SET	
MC4L-6S	8.5F X 15CM MULTI-CATH, QUAD LUMEN	10/BOX
MC4L-8S	8.5F X 20CM MULTI-CATH, QUAD LUMEN	10/BOX

- BASIC SET CONTAINS:**
- (1) Catheter w/ Clamps
 - (1) Guidewire w/ Advancer
 - (1) Introducer Needle
 - (1 | 2 | 3 | 4) Injection Cap(s)
 - (1) 5cc Syringe
 - (1) Vessel Dilator
 - (1) Movable Wing