



**Polyurethane**

Material

**Split Tip**

Tip Design

**15F**

French Size

**Straight**

Configuration

**Antegrade**

Insertion Type



**I-BEAM INTERNAL LUMEN DESIGN**  
Provides superior kink resistance  
and high flow rates.

BASIC SET		
JDDF1524	15F X 24CM ESCHOLON	5/BOX
JDDF1528	15F X 28CM ESCHOLON	5/BOX
JDDF1532	15F X 32CM ESCHOLON	5/BOX
JDDF1536	15F X 36CM ESCHOLON	5/BOX
JDDF1540	15F X 40CM ESCHOLON	5/BOX

**BASIC SET CONTAINS:**  
 (1) Catheter w/ Stylet  
 (1) Introducer Needle  
 (1) Guidewire  
 (2) Dilators  
 (1) Scalpel  
 (1) Tunneler  
 (1) Peelable Sheath Introducer  
 (1) Adhesive Wound Dressing  
 (2) Injection Caps



**SPLIT TIP DESIGN**  
Improves patency, flow  
and lowers recirculation.

**TWO-PIECE EXTENSION**  
Design enables the easy  
replacement of catheter  
extensions.

