

SPLIT CATH® RG


LONG TERM



Polyurethane	Split Tip	14F	Straight	Antegrade
Material	Tip Design	French Size	Configuration	Insertion Type



DOUBLE "D" INTERNAL LUMEN DESIGN



IMPLANTABLE CUFF

SPLITTABLE TIP LUMENS
Two free-floating lumens improve patency, flow and lowers recirculation

EASY SNAP HUB
Gives a secure and rapid attachment to the catheter after back-tunneling.

	BASIC SET	
MR771401	14F X 24CM SPLIT CATH® RG	5/BOX
MR771402	14F X 28CM SPLIT CATH® RG	5/BOX
MR771403	14F X 32CM SPLIT CATH® RG	5/BOX
MR771404	14F X 36CM SPLIT CATH® RG	5/BOX
MR771405	14F X 40CM SPLIT CATH® RG	5/BOX

BASIC SET CONTAINS:

- (1) Catheter w/ Stylet
- (1) Extension
- (1) Hub Clamp Assembly
- (2) Compression Rings
- (1) Lumen Clamp
- (1) Raulerson Ring Handled Tunneler
- (1) 5" Triball Tunneler
- (2) Dilators
- (1) Valved Peelable Introducer
- (1) Marked Guidewire
- (1) Scalpel
- (2) Injection Caps
- (1) Adhesive Wound Dressing
- (1) Patient ID Label