




DOUBLE "D" INTERNAL LUMEN DESIGN

Polyurethane	Split Tip	14F, 16F	Straight	Antegrade
Material	Tip Design	French Sizes	Configuration	Insertion Type



RETROGRADE INSERTION
Ensures accurate tip placement, improving dialysis adequacy.



REMOVABLE HUB
Aids in long-term patient comfort.



SPLIT TIP DESIGN
Improves patency, flow and lowers recirculation.

	BASIC SET	
SST24	14F X 24CM SPLIT STREAM®	5/BOX
SST28	14F X 28CM SPLIT STREAM®	5/BOX
SST32	14F X 32CM SPLIT STREAM®	5/BOX
SST36	14F X 36CM SPLIT STREAM®	5/BOX
SST40	14F X 40CM SPLIT STREAM®	5/BOX
SST2416	16F X 24CM SPLIT STREAM®	5/BOX
SST2816	16F X 28CM SPLIT STREAM®	5/BOX
SST3216	16F X 32CM SPLIT STREAM®	5/BOX
SST3616	16F X 36CM SPLIT STREAM®	5/BOX
SST4016	16F X 40CM SPLIT STREAM®	5/BOX

BASIC SET CONTAINS:

- (1) Catheter
- (1) Guidewire
- (2) Dilators
- (1) Arterial Extension Set (Red)
- (1) Venous Extension Set (Blue)
- (1) Scalpel
- (1) Valved Peelable Introducer
- (1) Adhesive Wound Dressing
- (4) Clamps
- (1) Removable Hub
- (1) Tunneler
- (1) Lumen Clamp
- (2) Injection Caps
- (1) Introducer Needle
- (1) Catheter Securement Device
- (1) Patient I.D. Card