




INDEPENDENT SINGLE ROUND INTERNAL LUMEN DESIGN

<b>Polyurethane</b> Material	<b>Single Round</b> Tip Design	<b>6.5F, 10F</b> French Sizes	<b>Straight</b> Configuration	<b>Retrograde</b> Insertion Type
---------------------------------	-----------------------------------	----------------------------------	----------------------------------	-------------------------------------

RETROGRADE INSERTION Ensures accurate tip placement, improving dialysis adequacy.



FOOTBALL CUFF Secures device for long-term use.

TWO INDEPENDENT, ROUND LUMENS Provides unparalleled patency and high flow characteristics.

BASIC SET		
PBF	6.5F X 29CM TESIO®, CUFFS @ 12CM AND 15CM	5/BOX
BFS-6	10F X 52CM TESIO®, CUFFS @ 18CM AND 21CM	5/BOX
BFR-6	10F X 52CM TESIO®, CUFFS @ 22CM AND 25CM	5/BOX
BFL-6	10F X 52CM TESIO®, CUFFS @ 27CM AND 30CM	5/BOX
BFR1035KD-A	10F X 40CM TESIO®, CUFFS @ 18CM AND 22CM	5/BOX
BFR1070KD-A	10F X 70CM TESIO®, CUFFS @ 46CM AND 50CM	5/BOX

BASIC SET WITH CUTTING TROCAR		
BFS-6CTT	10F X 52CM TESIO®, CUFFS @ 18CM AND 21CM	5/BOX
BFR-6CTT	10F X 52CM TESIO®, CUFFS @ 22CM AND 25CM	5/BOX
BFL-6CTT	10F X 52CM TESIO®, CUFFS @ 27CM AND 30CM	5/BOX

SINGLE SET		
BFSS	10F X 52CM TESIO®, CUFF @ 21CM	5/BOX
BFRS	10F X 52CM TESIO®, CUFF @ 25CM	5/BOX
BFLS	10F X 52CM TESIO®, CUFF @ 30CM	5/BOX

ACCESSORIES		
BTT	TESIO® BLUNT TUNNELER	5/BOX
CTT	10F TESIO® CUTTING TROCAR	5/BOX
FTT10	TESIO® FOOTBALL TUNNELER	5/BOX
PCTT	6.5F TESIO® CUTTING TROCAR	5/BOX
RTT	RAULERSON TUNNELER	5/BOX
BF11A	10F ARTERIAL EXTENSION SET FOR TESIO® CATHETER	10/BOX
BF11V	10F VENOUS EXTENSION SET FOR TESIO® CATHETER	10/BOX
BF11AVS	10F ARTERIAL AND VENOUS EXTENSION SET FOR TESIO® CATHETER	10/BOX
PBF11A	6.5F ARTERIAL EXTENSION SET FOR PEDIATRIC TESIO® CATHETER	10/BOX
PBF11V	6.5F VERNOUS EXTENSION SET FOR PEDIATRIC TESIO® CATHETER	10/BOX

BASIC SET CONTAINS:

- (2) Catheters
- (2) Introducer Needles
- (2) Guidewires
- (1) Scalpel
- (2) Tunnelers
- (2) Football Tunnelers
- (1) Coaxial Introducer Sheath
- (2) Valved Peelable Introducers
- (1) Arterial Extension Set (Red)
- (1) Venous Extension Set (Blue)
- (1) Adhesive Wound Dressing
- (2) Clamps
- (2) Injection Caps
- (1) Patient I.D. Card

BASIC SET WITH CUTTING TROCAR CONTAINS:

- (2) Catheters
- (2) Introducer Needles
- (2) Guidewires
- (1) Scalpel
- (2) Cutting Trocars
- (2) Football Tunnelers
- (1) Coaxial Introducer Sheath
- (2) Valved Peelable Introducers
- (1) Arterial Extension Set (Red)
- (1) Venous Extension Set (Blue)
- (1) Adhesive Wound Dressing
- (2) Clamps
- (2) Injection Caps
- (1) Patient I.D. Card

SINGLE SET CONTAINS:

- (1) Catheter
- (1) Introducer Needle
- (1) Guidewire
- (1) Scalpel
- (1) Tunneler
- (1) Football Tunneler
- (1) Peelable Sheath Introducer
- (1) Venous Extension
- (1) Adhesive Wound Dressing
- (1) Clamp
- (1) Injection Cap
- (1) Patient I.D. Card