



Polyurethane Material	Step Tip Tip Design	15.5F French Size	Straight, Pre-Curved Configurations	Antegrade Insertion Type
---------------------------------	-------------------------------	-----------------------------	---	------------------------------------



DOUBLE "D" INTERNAL LUMEN DESIGN



TAPERED, FAST TRACK STYLET
Provides a less traumatic and faster insertion.

90 DEGREE ARTERIAL TIP
Reduces risk of arterial insufficiency and enables full heparin lock.

90° Degree

15.5F LUMEN
Enables high flow rates.

BASIC SET WITH SIDHOLES		
THD155024	15.5F X 24CM TITAN™, STRAIGHT WITH SIDHOLES	5/BOX
THD155028	15.5F X 28CM TITAN™, STRAIGHT WITH SIDHOLES	5/BOX
THD155032	15.5F X 32CM TITAN™, STRAIGHT WITH SIDHOLES	5/BOX
THD155036	15.5F X 36CM TITAN™, STRAIGHT WITH SIDHOLES	5/BOX
THD155040	15.5F X 40CM TITAN™, STRAIGHT WITH SIDHOLES	5/BOX
THD155055	15.5F X 55CM TITAN™, STRAIGHT WITH SIDHOLES	5/BOX
THD155224	15.5F X 24CM TITAN™, PRE-CURVED WITH SIDHOLES	5/BOX
THD155228	15.5F X 28CM TITAN™, PRE-CURVED WITH SIDHOLES	5/BOX
THD155232	15.5F X 32CM TITAN™, PRE-CURVED WITH SIDHOLES	5/BOX
THD155236	15.5F X 36CM TITAN™, PRE-CURVED WITH SIDHOLES	5/BOX

BASIC SET CONTAINS:

- (1) Catheter w/ Stylet
- (1) Introducer Needle
- (1) Guidewire
- (2) Dilators
- (1) Scalpel
- (1) Tunneler
- (1) Valved Peelable Introducer
- (1) Adhesive Wound Dressing
- (2) Injection Caps
- (1) Patient I.D. Card

BASIC SET WITHOUT SIDHOLES		
THD155424	15.5F X 24CM TITAN™, STRAIGHT WITHOUT SIDHOLES	5/BOX
THD155428	15.5F X 28CM TITAN™, STRAIGHT WITHOUT SIDHOLES	5/BOX
THD155432	15.5F X 32CM TITAN™, STRAIGHT WITHOUT SIDHOLES	5/BOX
THD155436	15.5F X 36CM TITAN™, STRAIGHT WITHOUT SIDHOLES	5/BOX
THD155440	15.5F X 40CM TITAN™, STRAIGHT WITHOUT SIDHOLES	5/BOX
THD155455	15.5F X 55CM TITAN™, STRAIGHT WITHOUT SIDHOLES	5/BOX
THD155624	15.5F X 24CM TITAN™, PRE-CURVED WITHOUT SIDHOLES	5/BOX
THD155628	15.5F X 28CM TITAN™, PRE-CURVED WITHOUT SIDHOLES	5/BOX
THD155632	15.5F X 32CM TITAN™, PRE-CURVED WITHOUT SIDHOLES	5/BOX
THD155636	15.5F X 36CM TITAN™, PRE-CURVED WITHOUT SIDHOLES	5/BOX