



**Single,  
Double**

Lumens

**60cm**

Length

**3F, 4F, 5F**

French Sizes



### CVP MONITORING

All valved PICCs are indicated for Central Venous Pressure Monitoring.

### SALINE ONLY FLUSH

Limits the potential complications & costs of heparin use.

### CLAMP FREE EXTENSIONS

Eliminates maintenance challenges associated with needless injection caps.

CATHETER ONLY		
MR82013100	3F X 60CM VASCU-PICC® VALVED, SINGLE LUMEN	10/BOX
MR82014100	4F X 60CM VASCU-PICC® VALVED, SINGLE LUMEN	10/BOX
MR82014200	4F X 60CM VASCU-PICC® VALVED, DOUBLE LUMEN	10/BOX
MR82015200	5F X 60CM VASCU-PICC® VALVED, DOUBLE LUMEN	10/BOX

BASIC SET		
MR82013101	3F X 60CM VASCU-PICC VALVED, SINGLE LUMEN	5/BOX
MR82014101	4F X 60CM VASCU-PICC VALVED, SINGLE LUMEN	5/BOX
MR82014201	4F X 60CM VASCU-PICC VALVED, DOUBLE LUMEN	5/BOX
MR82015201	5F X 60CM VASCU-PICC VALVED, DOUBLE LUMEN	5/BOX

NURSING KIT		
MR82013105	3F X 60CM VASCU-PICC® VALVED, SINGLE LUMEN	5/BOX
MR82014105	4F X 60CM VASCU-PICC® VALVED, SINGLE LUMEN	5/BOX
MR82014205	4F X 60CM VASCU-PICC® VALVED, DOUBLE LUMEN	5/BOX
MR82015205	5F X 60CM VASCU-PICC® VALVED, DOUBLE LUMEN	5/BOX

FULL NURSING KIT		
MR82013106	3F X 60CM SINGLE VASCU-PICC® FULL NURSING	5/BOX
MR82014106	4F X 60CM SINGLE VASCU-PICC® FULL NURSING	5/BOX
MR82014206	4F X 60CM DUAL VASCU-PICC® FULL NURSING	5/BOX
MR82015206	5F X 60CM DUAL VASCU-PICC® FULL NURSING	5/BOX

MAXIMUM BARRIER		
MR82014106MB	4F X 60CM SINGLE VASCU-PICC® MAX BARRIER	3/BOX
MR82015206MB	5F X 60CM DUAL VASCU-PICC® MAX BARRIER	3/BOX

#### CATH ONLY CONTAINS:

- (1) PICC
- (1|2) Needle Free Valve(s)
- (1) Patient Information Pack
- (1) Patient Chart Sticker
- (1) Catheter Securement Device
- (1) Twisted Wire Stylet

#### BASIC SET CONTAINS:

- (1) PICC
- (1) Peelable Sheath Introducer
- (1) Mini Scalpel
- (1) Guidewire
- (1) 10cc Syringe
- (1) Tape Measure
- (1) Introducer Needle
- (1|2) Needle Free Valve(s)
- (1) Catheter Securement Device
- (1) Patient Information Pack
- (1) Patient Chart Sticker

#### NURSING KIT CONTAINS:

- (1) PICC
- (1) Peelable Sheath Introducer
- (1) Safety Scalpel
- (1) Guidewire
- (1) 10cc Syringe
- (1) Tape Measure
- (1) Safety Introducer Needle
- (1|2) Needle Free Valve(s)
- (1) Catheter Securement Device
- (1) Latex-Free Tourniquet
- (1) Patient Information Pack
- (1) Patient Chart Sticker
- (1) Twisted Wire Stylet

#### FULL NURSING KIT CONTAINS:

- (1) PICC
- (1) Peelable Sheath Introducer
- (1) Twisted Wire Stylet
- (1) Safety Scalpel
- (1) Guidewire
- (1) 5cc Syringe
- (2) Tape Measures
- (1) Surgical Tape
- (1) Safety Introducer Needle
- (1|2) Needle Free Valve(s)
- (1) Catheter Securement Device
- (1) Latex-Free Tourniquet
- (3) Pre-Filled Saline Syringes
- (1) Chlorhexidine
- (1) Tinted Chlorhexidine
- (1) 5cc Lidocaine
- (1) Adhesive Wound Dressing
- (1) Needle Stick Pad
- (10) 4" x 4" Gauze
- (2) Surgical Gloves

#### MAXIMUM BARRIER

- FULL NURSING KIT CONTAINS:
- (1) PICC
- (1) Peelable Sheath Introducer
- (1) Safety Scalpel
- (1) Guidewire
- (1) 5cc Syringe
- (2) Tape Measures
- (1) Surgical Tape
- (1) Safety Introducer Needle
- (1|2) Needle Free Valve(s)
- (1) Catheter Securement Device
- (1) Latex-Free Tourniquet
- (3) Pre-Filled Saline Syringes
- (1) Tinted Chlorhexidine
- (1) Chlorhexidine
- (1) 5cc Lidocaine
- (1) Adhesive Wound Dressing
- (1) Needle Stick Pad
- (10) 4" x 4" Gauze
- (2) Surgical Gloves

- (1) Filter Straw
- (1) Scissors
- (1) Steri-strips
- (2) Face Masks
- (8) 2" x 2" Gauze
- (1) Towel
- (1) Fenestrated Drape
- (1) Drape
- (1) Hypodermic Safety Needle
- (1) Lidocaine/Saline Stickers
- (1) Full Body Drape
- (1) XL Gown
- (1) Surgical Face Mask w/ Eye Shield
- (1) Bouffant Cap
- (1) Saline Wipes
- (1) Introducer Needle
- (1) Patient Information Pack
- (1) Patient Chart Sticker
- (1) Twisted Wire Stylet