

Single, Double, Triple Lumens	60cm Length	3F, 4F, 5F French Sizes
---	-----------------------	-----------------------------------

REVERSE TAPER LUMEN DESIGN
Reduces post insertion bleeding and provides improved catheter strength against catheter kinking.

EASY-TO-READ DEPTH MARKS
Easily and accurately trim to the desired length.

EASY GRIP LUERS
Crack resistant luers combine comfort and strength.



CATHETER ONLY		
MR17013100	3F X 60CM VASCU-PICC®, SINGLE LUMEN	10/BOX
MR17014100	4F X 60CM VASCU-PICC®, SINGLE LUMEN	10/BOX
MR17014200	4F X 60CM VASCU-PICC®, DOUBLE LUMEN	10/BOX
MR17015200	5F X 60CM VASCU-PICC®, DOUBLE LUMEN	10/BOX
MR17016200	6F X 60CM VASCU-PICC®, DOUBLE LUMEN	10/BOX

BASIC SET		
MR17013101	3F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX
MR17014101	4F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX
MR17015101	5F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX
MR17014201	4F X 60CM VASCU-PICC®, DOUBLE LUMEN	5/BOX
MR17015201	5F X 60CM VASCU-PICC®, DOUBLE LUMEN	5/BOX
MR17016201	6F X 60CM VASCU-PICC®, DOUBLE LUMEN	5/BOX
MR17015301	5F X 60CM VASCU-PICC®, TRIPLE LUMEN	5/BOX

LONG WIRE SET		
MR17014102	4F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX
MR17015102	5F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX
MR17015202	5F X 60CM VASCU-PICC®, DOUBLE LUMEN	5/BOX
MR17016202	6F X 60CM VASCU-PICC®, DOUBLE LUMEN	5/BOX
MR17015302	5F X 60CM VASCU-PICC®, TRIPLE LUMEN	5/BOX

NURSING KIT		
MR17013105	3F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX
MR17014105	4F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX
MR17015105	5F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX
MR17014205	4F X 60CM VASCU-PICC®, DOUBLE LUMEN	5/BOX
MR17015205	5F X 60CM VASCU-PICC®, DOUBLE LUMEN	5/BOX
MR17015305	5F X 60CM VASCU-PICC®, TRIPLE LUMEN	5/BOX

FULL NURSING KIT		
MR17013106	3F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX
MR17014106	4F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX
MR17015106	5F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX
MR17014206	4F X 60CM VASCU-PICC®, DOUBLE LUMEN	5/BOX
MR17015206	5F X 60CM VASCU-PICC®, DOUBLE LUMEN	5/BOX
MR17015306	5F X 60CM VASCU-PICC®, TRIPLE LUMEN	5/BOX

MAXIMUM BARRIER KIT		
MR17014106MB	4F X 60CM VASCU-PICC®, SINGLE LUMEN	3/BOX
MR17015206MB	5F X 60CM VASCU-PICC®, DOUBLE LUMEN	3/BOX
MR17015306MB	5F X 60CM VASCU-PICC®, TRIPLE LUMEN	3/BOX

FLUSHABLE STYLET SET		
MR17014111	4F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX
MR17015111	5F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX

CATHETER ONLY CONTAINS:
 (1) PICC w/ Stylet
 (1) 2 Needle Free Valve(s)
 (1) Patient Information Pack
 (1) Patient Chart Sticker
 (1) Catheter Securement Device

BASIC SET CONTAINS:
 (1) PICC w/ Stylet
 (1) Peelable Sheath Introducer
 (1) Mini Scalpel
 (1) Guidewire
 (1) 10cc Syringe
 (1) Tape Measure
 (1) Introducer Needle
 (1) 2|3 Needle Free Valve(s)
 (1) Catheter Securement Device
 (1) Patient Information Pack
 (1) Patient Chart Sticker

LONG WIRE SET CONTAINS:
 (1) PICC
 (1) Peelable Sheath Introducer
 (1) Mini Scalpel
 (1) Guidewire
 (1) 10cc Syringe
 (1) Tape Measure
 (1) Introducer Needle
 (1) 2|3 Needle Free Valve(s)
 (1) Catheter Securement Device
 (1) Patient Information Pack
 (1) Patient Chart Sticker

NURSING KIT CONTAINS:
 (1) PICC w/ Stylet
 (1) Peelable Sheath Introducer
 (1) Safety Scalpel
 (1) Guidewire
 (1) 10cc Syringe
 (1) Tape Measure
 (1) Safety Introducer Needle
 (1) 2|3 Needle Free Valve(s)
 (1) Catheter Securement Device
 (1) Patient Information Pack
 (1) Patient Chart Sticker

FULL NURSING KIT CONTAINS:
 (1) PICC w/ Stylet
 (1) Peelable Sheath Introducer
 (1) Safety Scalpel
 (1) Guidewire
 (1) 5cc Syringe
 (2) Tape Measures
 (1) Surgical Tape
 (1) Safety Introducer Needle
 (1) 2|3 Needle Free Valve(s)
 (1) Catheter Securement Device
 (1) Latex-Free Tourniquet
 (3) Pre-Filled Saline Syringes
 (1) Chlorhexidine
 (1) Tinted Chlorhexidine
 (1) 5cc Lidocaine
contents continued...

(1) Adhesive Wound Dressing
 (1) Needle Stick Pad
 (10) 4" x 4" Gauze
 (3) Surgical Gloves
 (1) Filter Straw
 (1) Scissors
 (1) Steri-strips
 (2) Face Masks
 (8) 2" x 2" Gauze
 (1) I.V. Catheter
 (1) Towel
 (1) Fenestrated Drape
 (1) Drape
 (1) Safety Hypodermic Needle
 (1) Lidocaine/Saline Stickers
 (1) Forceps
 (1) Introducer Needle
 (1) Patient Information Pack
 (1) Patient Chart Sticker

MAXIMUM BARRIER KIT CONTAINS:
 (1) PICC w/ Stylet
 (1) Peelable Sheath Introducer
 (1) Safety Scalpel
 (1) Guidewire
 (1) 5cc Syringes
 (2) Tape Measures
 (1) Surgical Tape
 (1) Safety Introducer Needle
 (1) 2|3 Needle Free Valve(s)
 (1) Catheter Securement Device
 (1) Latex-Free Tourniquet
 (3) Pre-Filled Saline Syringes
 (1) Chlorhexidine
 (1) Tinted Chlorhexidine
 (1) 5cc Lidocaine
 (1) Adhesive Wound Dressing
 (1) Needle Stick Pad
 (10) 4" x 4" Gauze
 (3) Surgical Gloves
 (1) Filter Straw
 (1) Scissors
 (1) Steri-strips
 (2) Face Masks
 (8) 2" x 2" Gauze
 (1) I.V. Catheter
 (1) Towel
 (1) Fenestrated Drape
 (1) Drape
 (1) Safety Hypodermic Needle
 (1) Lidocaine/Saline Stickers
 (1) Full Body Drape
 (1) XL Gown
 (1) Surgical Face Mask w/ Eye Shield
 (1) Bouffant Cap
 (1) Introducer Needle
 (1) Saline Wipes
 (1) Forceps
 (1) Patient Information Pack
 (1) Patient Chart Sticker
 (1) Probe Cover

FLUSHABLE STYLET SET CONTAINS:
 *Contents same as BASIC SET