

PRO-FUSE®

FEATURES & BENEFITS

- Ports are approved for both CECT injections and infusion therapy.
- Allows CT injections for diagnostic imaging at up to 5cc/sec at 300psi.
- Raised rim provides clear tactile response during palpation to assist in septum location.
- Unique printing allows “CT” to be visible under X-ray.

Open Suture Holes

Catalog #	Pro-Fuse® CT Low Profile		
MRCTT80041	8.0F POLYURETHANE CATHETER	ATTACHABLE	5/Box
MRCTT66041	6.6F POLYURETHANE CATHETER	ATTACHABLE	5/Box
MRCTT80841	8.0F POLYURETHANE CATHETER	PRE-ATTACHED	5/Box
MRCTT66841	6.6F POLYURETHANE CATHETER	PRE-ATTACHED	5/Box

Silicone Filled Suture Holes

Catalog #	Pro-Fuse® CT		
MRCTT96001	9.6F SILICONE CATHETER	ATTACHABLE	5/Box
MRCTT80001	8.0F POLYURETHANE CATHETER	ATTACHABLE	5/Box
MRCTT96801	9.6F SILICONE CATHETER	PRE-ATTACHED	5/Box
MRCTT80801	8.0F POLYURETHANE CATHETER	PRE-ATTACHED	5/Box
	Pro-Fuse® CT Low Profile		
MRCTT8004S	8.0F POLYURETHANE CATHETER	ATTACHABLE	5/Box
MRCTT6604S	6.6F POLYURETHANE CATHETER	ATTACHABLE	5/Box
MRCTT8084S	8.0F POLYURETHANE CATHETER	PRE-ATTACHED	5/Box
MRCTT6684S	6.6F POLYURETHANE CATHETER	PRE-ATTACHED	5/Box

SET CONTAINS:

- (2) Catheter Locks (1) Scalpel (1) Introducer Needle (1) Vein Pick (1) Marked Guidewire (1) 10cc Syringe (1) Valved Peelable Introducer (1) Tuner (2) 22 Ga. Huber Needle (1 straight, 1 right angle) (1) Blunt Tip Needle (1) Valved Peelable Introducer IFU (1) Port IFU (1) Patient Information Pack (1) Patient Chart Sticker

SET CONTAINS:

- (2) Catheter Locks (1) Scalpel (1) Introducer Needle (1) Vein Pick (1) Marked Guidewire (1) 10cc Syringe (1) Valved Peelable Introducer (1) Tuner (2) 22 Ga. Huber Needle (1 straight, 1 right angle) (1) Blunt Tip Needle (1) Valved Peelable Introducer IFU (1) Port IFU (1) Patient Information Pack (1) Patient Chart Sticker



MR Conditional – 3 Tesla (artifacts may present imaging problems if MRI area of interest is on or near area where device is located)

