


NEW!

THE NEW, RETROGRADE Split Cath RG[®]

Using OTW Technique

The “**Over-The-Wire**” (OTW) approach eliminates the need for large peel-away-sheath introduction, which is especially helpful in left-sided internal jugular approaches. Experience the benefits of retrograde tip placement and the Split Cath[®] advantage.



- **Retrograde Approach** – Catheter positioning and tunneling are simplified.
- **Raulerson Tunneler** – Simplifies retrograde tunneling.
- **Benson Wire** – The rigid coil and Teflon coated surface greatly facilitates a stable Over-The-Wire approach.

Easy Snap Hub

Gives a secure and rapid attachment to the catheter after back-tunneling.

Splittable Tip Lumens

Increases patency to the level of single lumen catheters and Tesio catheters.

Multiple tunneling options providing versatility

Implantable Cuff

CATALOG CODE	DESCRIPTION	UOM
ASPC24OTW	14F X 24CM SPLIT CATH RG [®]	5/BOX
ASPC28OTW	14F X 28CM SPLIT CATH RG [®]	5/BOX
ASPC32OTW	14F X 32CM SPLIT CATH RG [®]	5/BOX
ASPC36OTW	14F X 36CM SPLIT CATH RG [®]	5/BOX
ASPC40OTW	14F X 40CM SPLIT CATH RG [®]	5/BOX

Set Contains:

- (1) Catheter
- (1) Venous/Arterial Extension Assembly
- (1) Attachable Hub
- (2) Injection Caps
- (1) Tunneling Tool
- (1) Guidewire
- (2) Vessel Dilators
- (1) Introducer Needle
- (1) Scalpel



1499 Delp Dr. Harleysville, PA 19438
P: 215-256-4201 F: 215-256-1787

www.medcompnet.com

PN2573 Rev. C 4/17