



**Plastic**  
Port Material

**Mini,  
Low Profile**  
Configurations

**Polyurethane**  
Catheter Material

**5cc/sec  
@ 325psi**  
Injection Rating

**Attachable,  
Pre-Attached**  
Catheter  
Configurations


**5F**  
Catheter Size



**REDUCED WIDTH AND  
SMALLER PROFILE**  
Smallest port for mini-  
mal port pocket and  
smaller incision size.



**CHEST OR ARM PLACEMENT**  
Provides alternative options for  
patients.



**PATIENT-CENTRIC DESIGN**  
Smooth contours and a tapered  
tail help minimize incision size  
while improving aesthetics and  
patient comfort.



PEDIATRIC KITS		
MRCTI50041P	5F DIGNITY® LP SHAPED ATTACHABLE	5/BOX
MRCTI50841P	5F DIGNITY® LP SHAPED PRE-ATTACHED	5/BOX
MRCTI50041MP	5F DIGNITY® MINI SHAPED ATTACHABLE	5/BOX
MRCTI50841MP	5F DIGNITY® MINI SHAPED PRE-ATTACHED	5/BOX
PEDIATRIC KITS W/ MICRO KIT		
MICTI50041P	5F DIGNITY® LP SHAPED ATTACHABLE	5/BOX
MICTI50841P	5F DIGNITY® LP SHAPED PRE-ATTACHED	5/BOX
MICTI50041MP	5F DIGNITY® MINI SHAPED ATTACHABLE	5/BOX
MICTI50841MP	5F DIGNITY® MINI SHAPED PRE-ATTACHED	5/BOX

**SET CONTAINS:**

- (1) Port
- (1) Catheter
- (2) Catheter Locks
- (1) Scalpel
- (1) Introducer Needle
- (1) Vein Pick
- (1) Marked Guidewire
- (1) 10cc Syringe
- (1) 3cc Syringe
- (1) Valved Peelable Introducer
- (1) Tunneler
- (2) 22 Ga. Huber Needles  
(1 straight, 1 right angle)
- (1) Blunt Tip Needle
- (1) Patient Information Pack
- (1) Patient Chart Sticker

**MICRO KIT CONTAINS:**

- (1) 5F Coaxial Micro-Stick®
- (1) Guidewire with  
Radiopaque Tip
- (1) 21ga Introducer Needle  
with Echogenic Tip