




Plastic Port Material	Mini Configuration	Polyurethane Catheter Material	5cc/sec @ 325psi Injection Rating	Attachable, Pre-Attached Catheter Configurations	5F, 6.6F, 8F French Sizes
---------------------------------	------------------------------	--	---	--	-------------------------------------



REDUCED WIDTH AND SMALLER PROFILE
Smallest port for minimal port pocket.



CHEST OR ARM PLACEMENT
Provides alternative options for patients.



PATIENT-CENTRIC DESIGN
Smooth contours and a tapered tail help minimize incision size while improving aesthetics and patient comfort.



DIGNITY® MINI	OPEN SUTURE HOLES		SILICONE FILLED SUTURE HOLES		QTY
		WITH MICRO INTRODUCER		WITH MICRO INTRODUCER	
5F ATTACHABLE	MRDP50ASN	MIDP50ASN	MRDP50AXS	MIDP50AXS	5/BOX
5F PRE-ATTACHED	MRDP50PSN	MIDP50PSN	MRDP50PSS	MIDP50PSS	5/BOX
6.6F ATTACHABLE	MRDP66ASN	MIDP66ASN	MRDP66AXS	MIDP66AXS	5/BOX
6.6F PRE-ATTACHED	MRDP66PSN	MIDP66PSN	MRDP66PSS	MIDP66PSS	5/BOX
8F ATTACHABLE	MRDP80ASN	MIDP80ASN	MRDP80AXS	MIDP80AXS	5/BOX
8F PRE-ATTACHED	MRDP80PSN	MIDP80PSN	MRDP80PSS	MIDP80PSS	5/BOX
DIGNITY MINI PEDIATRIC					
5F ATTACHABLE	MRCTI50041MP	MICTI50041MP	-	-	5/BOX
5F PRE-ATTACHED	MRCTI50841MP	MICTI50841MP	-	-	5/BOX

SET CONTAINS:

- (1) Port
- (1) Catheter
- (2) Catheter Locks
- (1) Scalpel
- (1) Introducer Needle
- (1) Vein Pick
- (1) Marked Guidewire
- (1) 10cc Syringe
- (1) Valved Peelable Introducer
- (1) Tunneler
- (2) 22ga Huber Needles
(1 Straight, 1 Right Angle)
- (1) Blunt Tip Needle
- (1) Patient Information Pack
- (1) Patient Chart Sticker

MICRO KIT CONTAINS:

- (1) 5F Coaxial Micro-Stick®
- (1) Guidewire with Radiopaque Tip
- (1) 21ga Introducer Needle with Echogenic Tip



MR Conditional – 3 Tesla (artifacts may present imaging problems if MRI area of interest is on or near area where device is located)