

Titanium Port Material	Mini, Low Profile, Mid Sized Configurations	Polyurethane, Silicone Catheter Materials	5cc/sec @ 325psi Injection Rating	Attachable, Pre-Attached Catheter Configurations	5F, 6.6F, 8F, 9.6F French Sizes
----------------------------------	---	---	---	--	---



TITANIUM BODY
Durable titanium material provides a more robust port.



CONCAVE SIDES
Allow for controlled, accurate port assembly and placement.



PATIENT-CENTRIC DESIGN
Smooth contours and a tapered tail help minimize incision size while improving aesthetics and patient comfort.

DIGNITY® MINI	OPEN SUTURE HOLES		SILICONE FILLED SUTURE HOLES		QTY
		WITH MICRO INTRODUCER		WITH MICRO INTRODUCER	
5F ATTACHABLE	MRDT50ASN	MIDT50ASN	MRDT50AXS	MIDT50AXS	5/BOX
5F PRE-ATTACHED	MRDT50PSN	MIDT50PSN	MRDT50PSS	MIDT50PSS	5/BOX
6.6F ATTACHABLE	MRDT66ASN	MIDT66ASN	MRDT66AXS	MIDT66AXS	5/BOX
6.6F PRE-ATTACHED	MRDT66PSN	MIDT66PSN	MRDT66PSS	MIDT66PSS	5/BOX
8F ATTACHABLE	MRDT80ASN	MIDT80ASN	MRDT80AXS	MIDT80AXS	5/BOX
8F PRE-ATTACHED	MRDT80PSN	MIDT80PSN	MRDT80PSS	MIDT80PSS	5/BOX
DIGNITY® LOW PROFILE					
5F ATTACHABLE	MRDT50ALN	MIDT50ALN	MRDT50ALS	MIDT50ALS	5/BOX
5F PRE-ATTACHED	MRDT50PLN	MIDT50PLN	MRDT50PLS	MIDT50PLS	5/BOX
6.6F ATTACHABLE	MRDT66ALN	MIDT66ALN	MRDT66ALS	MIDT66ALS	5/BOX
6.6F PRE-ATTACHED	MRDT66PLN	MIDT66PLN	MRDT66PLS	MIDT66PLS	5/BOX
8F ATTACHABLE	MRDT80ALN	MIDT80ALN	MRDT80ALS	MIDT80ALS	5/BOX
8F PRE-ATTACHED	MRDT80PLN	MIDT80PLN	MRDT80PLS	MIDT80PLS	5/BOX
DIGNITY® MID SIZED					
5F ATTACHABLE	MRDT50AMN	MIDT50AMN	MRDT50AMS	MIDT50AMS	5/BOX
5F PRE-ATTACHED	MRDT50PMN	MIDT50PMN	MRDT50PMS	MIDT50PMS	5/BOX
6.6F ATTACHABLE	MRDT66AMN	MIDT66AMN	MRDT66AMS	MIDT66AMS	5/BOX
6.6F PRE-ATTACHED	MRDT66PMN	MIDT66PMN	MRDT66PMS	MIDT66PMS	5/BOX
8F ATTACHABLE	MRDT80AMN	MIDT80AMN	MRDT80AMS	MIDT80AMS	5/BOX
8F PRE-ATTACHED	MRDT80PMN	MIDT80PMN	MRDT80PMS	MIDT80PMS	5/BOX
9.6F ATTACHABLE	MRDT96AMN	MIDT96AMN	MRDT96AMS	MIDT96AMS	5/BOX
9.6F PRE-ATTACHED	MRDT96PMN	MIDT96PMN	MRDT96PMS	MIDT96PMS	5/BOX

SETS CONTAIN:

- (1) Port
- (1) Catheter
- (2) Catheter Locks
- (1) Scalpel
- (1) Introducer Needle
- (1) Vein Pick
- (1) Marked Guidewire with Advancer
- (2) 10cc Syringes
- (1) Valved Peelable Introducer
- (1) Tunneler
- (2) 22ga Huber Needles (1-Straight, 1-Right Angle)
- (1) Blunt Tip Needle
- (1) Patient Information Pack
- (1) Patient Chart Sticker

SETS WITH MICRO-STICK® CONTAIN:

- (1) Port
- (1) Catheter
- (2) Catheter Locks
- (1) Scalpel
- (1) Vein Pick
- (1) Marked Guidewire with Advancer
- (2) 10cc Syringes
- (1) Valved Peelable Introducer
- (1) Tunneler
- (2) 22ga Huber Needles (1-Straight, 1-Right Angle)
- (1) Blunt Tip Needle
- (1) 21ga Introducer Needle
- (1) .018" X 45cm Nitinol Wire
- (1) 5F Coaxial Dilator Assembly
- (1) Patient Information Pack
- (1) Patient Chart Sticker

 Silicone Catheters



MR Conditional – 3 Tesla (artifacts may present imaging problems if MRI area of interest is on or near area where device is located)