




Polyurethane
Material

Tapered Tip
Tip Design

9F, 11.5F
French Sizes

**Straight,
Pre-Curved,
Raulerson IJ**

Configurations



TAPERED TIP
Reduces vessel trauma and provides an easier insertion.

COAXIAL INTERNAL LUMEN DESIGN
Round venous lumen provides greater hemodynamic capabilities.

RAULERSON IJ
Reduces infection and improves patient comfort during treatment.



CATHETER ONLY		
XTP94CT	9F X 12CM DUO-FLOW®, STRAIGHT	10/BOX
XTP96CT	9F X 15CM DUO-FLOW®, STRAIGHT	10/BOX
XTP98CT	9F X 20CM DUO-FLOW®, STRAIGHT	10/BOX
XTP94IJC	9F X 12CM DUO-FLOW®, PRE-CURVED	10/BOX
XTP96IJC	9F X 15CM DUO-FLOW®, PRE-CURVED	10/BOX
XTP98IJC	9F X 20CM DUO-FLOW®, PRE-CURVED	10/BOX

BASIC SET		
XTP94MT	9F X 12CM DUO-FLOW®, STRAIGHT	5/BOX
XTP96MT	9F X 15CM DUO-FLOW®, STRAIGHT	5/BOX
XTP98MT	9F X 20CM DUO-FLOW®, STRAIGHT	5/BOX
XTP94IJS	9F X 12CM DUO-FLOW®, PRE-CURVED	5/BOX
XTP96IJS	9F X 15CM DUO-FLOW®, PRE-CURVED	5/BOX
XTP98IJS	9F X 20CM DUO-FLOW®, PRE-CURVED	5/BOX

CATHETER ONLY CONTAINS:

- (1) Catheter
- (2) Injection Caps
- (1) Dilator
- (1) Patient I.D. Card

BASIC SET CONTAINS:

- (1) Catheter
- (1) Introducer Needle
- (1) Guidewire
- (1) Dilator
- (1) Scalpel
- (1) Adhesive Wound Dressing
- (2) Injection Caps
- (1) Patient I.D. Card

FULL TRAY CONTAINS:

- (1) Catheter
- (1) Introducer Needle
- (1) Guidewire
- (1) Dilator
- (1) Safety Scalpel
- (1) Adhesive Wound Dressing
- (2) Injection Caps
- (8) 4" x 4" Gauze
- (2) Chlorhexidine
- (2) 5cc Syringes
- (2) Pre-Filled Saline Syringes
- (1) 2-0 Silk Suture
- (3) Safety Hypodermic Needles
- (1) Hemostat
- (1) Fenestrated Drape
- (1) 5cc Lidocaine
- (1) Filter Straw
- (1) Surgical Gloves
- (1) Needle Stick Pad
- (1) Patient I.D. Card

US and Foreign Patents Pending.

CATHETER ONLY		
XTP114CT	11.5F X 12CM DUO-FLOW®, STRAIGHT	10/BOX
XTP116CT	11.5F X 15CM DUO-FLOW®, STRAIGHT	10/BOX
XTP118CT	11.5F X 20CM DUO-FLOW®, STRAIGHT	10/BOX
XTP119CT	11.5F X 24CM DUO-FLOW®, STRAIGHT	10/BOX
XTP114IJC	11.5F X 12CM DUO-FLOW®, RAULERSON IJ	10/BOX
XTP116IJC	11.5F X 15CM DUO-FLOW®, RAULERSON IJ	10/BOX
XTP118IJC	11.5F X 20CM DUO-FLOW®, RAULERSON IJ	10/BOX
XTP114IJC-2	11.5F X 12CM DUO-FLOW®, PRE-CURVED	10/BOX
XTP116IJC-2	11.5F X 15CM DUO-FLOW®, PRE-CURVED	10/BOX
XTP118IJC-2	11.5F X 20CM DUO-FLOW®, PRE-CURVED	10/BOX

BASIC SET		
XTP114MT	11.5F X 12CM DUO-FLOW®, STRAIGHT	5/BOX
XTP116MT	11.5F X 15CM DUO-FLOW®, STRAIGHT	5/BOX
XTP118MT	11.5F X 20CM DUO-FLOW®, STRAIGHT	5/BOX
XTP119MT	11.5F X 24CM DUO-FLOW®, STRAIGHT	5/BOX
XTP114IJS	11.5F X 12CM DUO-FLOW®, RAULERSON IJ	5/BOX
XTP116IJS	11.5F X 15CM DUO-FLOW®, RAULERSON IJ	5/BOX
XTP118IJS	11.5F X 20CM DUO-FLOW®, RAULERSON IJ	5/BOX
XTP114IJS-2	11.5F X 12CM DUO-FLOW®, PRE-CURVED	5/BOX
XTP116IJS-2	11.5F X 15CM DUO-FLOW®, PRE-CURVED	5/BOX
XTP118IJS-2	11.5F X 20CM DUO-FLOW®, PRE-CURVED	5/BOX

FULL TRAY		
MCDLT3500	11.5F X 12CM DUO-FLOW®, STRAIGHT	5/BOX
MCDLT4000	11.5F X 15CM DUO-FLOW®, STRAIGHT	5/BOX
MCDLT4400	11.5F X 20CM DUO-FLOW®, STRAIGHT	5/BOX
MCDLT4424#	11.5F X 24CM DUO-FLOW®, STRAIGHT	5/BOX
MCDLT114IJ	11.5F X 12CM DUO-FLOW®, RAULERSON IJ	5/BOX
MCDLT116IJ	11.5F X 15CM DUO-FLOW®, RAULERSON IJ	5/BOX
MCDLT118IJ	11.5F X 20CM DUO-FLOW®, RAULERSON IJ	5/BOX
MCDLT114IJ-2	11.5F X 12CM DUO-FLOW®, PRE-CURVED	5/BOX
MCDLT116IJ-2	11.5F X 15CM DUO-FLOW®, PRE-CURVED	5/BOX
MCDLT118IJ-2	11.5F X 20CM DUO-FLOW®, PRE-CURVED	5/BOX