



Silicone Material	Tapered Tip Tip Design	8F, 11.5F French Sizes	Straight Configuration
-----------------------------	----------------------------------	----------------------------------	----------------------------------



DOUBLE "D" INTERNAL LUMEN DESIGN



FAST TRACK STYLET
Permits easier insertion.



SILICONE MATERIAL
Provides biocompatibility, kink resistance and improved patient comfort.



TAPERED TIP
Reduces vessel trauma and provides an easier insertion.



CATHETER ONLY		
SL15C	11.5F X 15CM HEMO-CATH® ST	10/BOX
SL20C	11.5F X 20CM HEMO-CATH® ST	10/BOX
SL24C	11.5F X 24CM HEMO-CATH® ST	10/BOX

8F BASIC SET		
SL12P	8F X 12CM HEMO-CATH® ST	5/BOX

11.5F BASIC SET		
SL15	11.5F X 15CM HEMO-CATH® ST	5/BOX
SL20	11.5F X 20CM HEMO-CATH® ST	5/BOX
SL24	11.5F X 24CM HEMO-CATH® ST	5/BOX
SL27	11.5F X 27CM HEMO-CATH® ST	5/BOX

FULL TRAY		
MCDLTSL15	11.5F X 15CM HEMO-CATH® ST	5/BOX
MCDLTSL20	11.5F X 20CM HEMO-CATH® ST	5/BOX
MCDLTSL24	11.5F X 24CM HEMO-CATH® ST	5/BOX
MCDLTSL27	11.5F X 27CM HEMO-CATH® ST	5/BOX

CATHETER ONLY CONTAINS:

- (1) Catheter w/ Stylet
- (2) Injection Caps
- (1) Hemo-Cath® Hub Clip
- (1) Patient I.D. Card

8F BASIC SET CONTAINS:

- (1) Catheter
- (1) Introducer Needle
- (1) Guidewire
- (1) Dilator
- (1) 10cc Syringe
- (1) Adhesive Wound Dressing
- (2) Injection Caps
- (1) Hemo-Cath® Hub Clip
- (1) Patient I.D. Card

11.5F BASIC SET CONTAINS:

- (1) Catheter w/ Stylet
- (1) Introducer Needle
- (1) Guidewire
- (2) Dilators
- (1) Adhesive Wound Dressing
- (2) Injection Caps
- (1) Hemo-Cath® Hub Clip
- (1) Patient I.D. Card

FULL TRAY CONTAINS:

- (1) Catheter w/ Stylet
- (1) Safety Introducer Needle
- (1) Guidewire
- (1) Dilator
- (1) Safety Scalpel
- (1) Peelable Sheath Introducer
- (1) Adhesive Wound Dressing
- (2) Injection Caps
- (1) Hemo-Cath® Hub Clip
- (4) 4" x 4" Gauze
- (2) Chlorhexidine
- (2) 5cc Syringe
- (2) Pre-Filled Saline Syringes
- (1) 2-0 Silk Suture
- (3) Safety Hypodermic Needles
- (1) Hemostat
- (1) Fenestrated Drape
- (1) 5cc Lidocaine
- (1) Filter Straw
- (1) Surgical Gloves
- (1) Needle Stick Pad
- (1) Patient I.D. Card