


| | | | |
|-----------------------------|----------------------------------|---------------------------|----------------------------------|
| Silicone Material | Tapered Tip Tip Design | 14F French Size | Straight Configuration |
|-----------------------------|----------------------------------|---------------------------|----------------------------------|



DOUBLE "D" INTERNAL LUMEN DESIGN



14F LUMEN
Allows for improved flow rates.



TAPERED TIP
Reduces vessel trauma and provides an easier insertion.



SILICONE MATERIAL
Provides biocompatibility, kink resistance and improved patient comfort.

| | CATHETER ONLY | |
|----------|----------------|-------|
| MC061442 | 14F X 15CM SLX | 5/BOX |
| MC061443 | 14F X 20CM SLX | 5/BOX |
| MC061444 | 14F X 24CM SLX | 5/BOX |
| MC061445 | 14F X 30CM SLX | 5/BOX |

| | BASIC SET | |
|----------|----------------|-------|
| MC061402 | 14F X 15CM SLX | 5/BOX |
| MC061403 | 14F X 20CM SLX | 5/BOX |
| MC061404 | 14F X 24CM SLX | 5/BOX |
| MC061405 | 14F X 30CM SLX | 5/BOX |

| | FULL TRAY | |
|----------|----------------|-------|
| MC061422 | 14F X 15CM SLX | 5/BOX |
| MC061423 | 14F X 20CM SLX | 5/BOX |
| MC061424 | 14F X 24CM SLX | 5/BOX |
| MC061425 | 14F X 30CM SLX | 5/BOX |

| | MAXIMUM BARRIER | |
|------------|-----------------|-------|
| MC061422MB | 14F X 15CM SLX | 4/BOX |
| MC061423MB | 14F X 20CM SLX | 4/BOX |
| MC061424MB | 14F X 24CM SLX | 4/BOX |
| MC061425MB | 14F X 30CM SLX | 4/BOX |

CATHETER ONLY CONTAINS:

- (1) Catheter w/ Stylet
- (2) Injection Caps
- (1) Hemo-Cath® Hub Clip
- (1) Patient I.D. Card

BASIC SET CONTAINS:

- (1) Catheter w/ Stylet
- (1) Introducer Needle
- (1) Guidewire
- (2) Dilators
- (1) Scalpel
- (1) Peelable Sheath Introducer
- (1) Adhesive Wound Dressing
- (2) Injection Caps
- (1) Hemo-Cath® Hub Clip
- (1) 2-0 Suture
- (1) Patient I.D. Card

FULL TRAY CONTAINS:

- (1) Catheter w/ Stylet
- (1) Safety Introducer Needle
- (1) Guidewire
- (2) Dilators
- (1) Safety Scalpel
- (1) Peelable Sheath Introducer
- (1) Adhesive Wound Dressing
- (2) Injection Caps
- (1) Hemo-Cath® Hub Clip
- (8) 2" x 2" Gauze
- (2) Chlorhexidine
- (1) 5cc Syringe
- (2) Pre-Filled Saline Syringes
- (1) 2-0 Silk Suture
- (3) Safety Hypodermic Needles
- (1) Hemostat
- (1) Fenestrated Drape
- (1) 5cc Lidocaine
- (1) Filter Straw
- (1) Surgical Gloves
- (1) Needle Stick Pad
- (1) Patient I.D. Card

MAXIMUM BARRIER SET CONTAINS:

- (1) Catheter w/ Stylet
- (1) Safety Introducer Needles
- (1) Guidewire
- (1) Guidewire Advancer
- (2) Dilators
- (1) Safety Scalpel
- (1) Peelable Sheath Introducer
- (1) Adhesive Wound Dressing
- (2) Injection Caps
- (1) Hemo-Cath Hub Clip
- (10) 4" x 4" Gauze
- (8) 2" x 2" Gauze
- (1) Chlorhexidine
- (1) 5cc Syringe
- (3) Pre-Filled Saline Syringes
- (1) 2-0 Silk Suture
- (3) Safety Hypodermic Needles
- (1) Hemostat
- (1) Fenestrated Drape
- (1) 5cc Lidocaine
- (1) Filter Straw
- (3) Surgical Gloves
- (1) Needle Stick Pad
- (1) Patient I.D. Card
- (1) Absorbant Towel
- (1) Steri-Strips
- (2) Face Mask
- (1) Lidocaine/Saline Stickers
- (1) Tinted Chlorhexidine
- (1) Surgical Tape
- (1) XL Gown
- (1) Bouffant Cap
- (1) Surgical Mask

