



**REMOVABLE “Y” ADAPTER**  
Also includes pre-loaded stylet.

**ROTATING SUTURE WING**  
Provides secure external anchoring.

**FLEXI-SMOOTH MATERIAL**  
Softens in-situ, conforming to vessel anatomy.

	CATHETER ONLY	
MCY306PS	8F X 15CM SUBCLAVIAN CATHETER, WITH REMOVABLE “Y” HUB	10/BOX
MCY308PS	8F X 20CM SUBCLAVIAN CATHETER, WITH REMOVABLE “Y” HUB	10/BOX

**CATHETER ONLY CONTAINS:**  
 (1) Catheter w/ Stylet  
 (1) Removable “Y” Adapter  
 (1) Injection Cap  
 (1) Patient I.D. Card