




Polyurethane Material	Step Tip Tip Design	15.5F French Size	Straight Configuration
---------------------------------	-------------------------------	-----------------------------	----------------------------------



TRIPLE INTERNAL LUMEN DESIGN



	CATHETER ONLY	
MC013541	15.5F X 15CM T-3 CT, STRAIGHT EXTENSION	10/BOX
MC013542	15.5F X 20CM T-3 CT, STRAIGHT EXTENSION	10/BOX
MC013543	15.5F X 24CM T-3 CT, STRAIGHT EXTENSION	10/BOX
MC013544	15.5F X 28CM T-3 CT, STRAIGHT EXTENSION	10/BOX
MC013545	15.5F X 32CM T-3 CT, STRAIGHT EXTENSION	10/BOX

	BASIC SET	
MC013501	15.5F X 15CM T-3 CT, STRAIGHT EXTENSION	5/BOX
MC013502	15.5F X 20CM T-3 CT, STRAIGHT EXTENSION	5/BOX
MC013503	15.5F X 24CM T-3 CT, STRAIGHT EXTENSION	5/BOX
MC013504	15.5F X 28CM T-3 CT, STRAIGHT EXTENSION	5/BOX
MC013505	15.5F X 32CM T-3 CT, STRAIGHT EXTENSION	5/BOX

	FULL TRAY	
MC013521	15.5F X 15CM T-3 CT, STRAIGHT EXTENSION	5/BOX
MC013522	15.5F X 20CM T-3 CT, STRAIGHT EXTENSION	5/BOX
MC013523	15.5F X 24CM T-3 CT, STRAIGHT EXTENSION	5/BOX
MC013524	15.5F X 28CM T-3 CT, STRAIGHT EXTENSION	5/BOX
MC013525	15.5F X 32CM T-3 CT, STRAIGHT EXTENSION	5/BOX

	MAXIMUM BARRIER	
MC013521MB	15.5F X 15CM T-3 CT, STRAIGHT EXTENSION	4/BOX
MC013522MB	15.5F X 20CM T-3 CT, STRAIGHT EXTENSION	4/BOX
MC013523MB	15.5F X 24CM T-3 CT, STRAIGHT EXTENSION	4/BOX
MC013524MB	15.5F X 28CM T-3 CT, STRAIGHT EXTENSION	4/BOX
MC013525MB	15.5F X 32CM T-3 CT, STRAIGHT EXTENSION	4/BOX

CATHETER ONLY CONTAINS:

- (1) Catheter w/ Stylet
- (3) Dilators
- (3) Injection Caps

BASIC SET CONTAINS:

- (1) Catheter w/ Stylet
- (1) Introducer Needle
- (1) Guidewire
- (3) Dilators
- (1) Scalpel
- (1) Valved Peelable Introducer
- (1) Adhesive Wound Dressing
- (3) Injection Caps

FULL TRAY CONTAINS:

- (1) Catheter w/ Stylet
- (1) Introducer Needle
- (1) Guidewire
- (3) Dilators
- (1) Safety Scalpel
- (1) Valved Peelable Introducer
- (1) Adhesive Wound Dressing
- (3) Injection Caps
- (4) 4" x 4" Gauze
- (2) Chlorhexidine
- (1) 5cc Syringe
- (3) 10cc Pre-Filled Saline Syringes
- (1) 2-0 Silk Suture
- (3) Safety Hypodermic Needles
- (1) Hemostat
- (1) Fenestrated Drape
- (1) 5cc Lidocaine
- (1) Filter Straw
- (1) Surgical Gloves

MAXIMUM BARRIER CONTAINS:

- (1) Catheter
- (1) Introducer Needle
- (1) Guidewire
- (3) Dilators
- (1) Safety Scalpel
- (1) Valved Introducer
- (1) Adhesive Wound Dressing
- (3) Injection Caps
- (4) 4" x 4" Gauze
- (8) 2" x 2" Gauze
- (2) Chlorhexidine
- (1) 5cc Syringe
- (3) 10cc Pre-Filled Syringes
- (1) 2-0 Silk Suture
- (3) Safety Hypodermic Needles
- (1) Hemostat
- (1) Fenestrated Drape
- (1) 5cc Lidocaine
- (1) Filter Straw
- (3) Pairs of Surgical Gloves
- (1) XL Drape
- (1) XL Gown
- (1) Tinted Chlorhexidine
- (1) Surgical Tape
- (1) Ear Loop Face Mask
- (1) Steri-Strips
- (1) Lidocaine/Saline Stickers
- (1) Bouffant Cap
- (1) Surgical Facemask