



TRIPLE INTERNAL LUMEN DESIGN

Polyurethane
Material

Tapered Tip
Tip Design

11.5F
French Size

**Straight,
Pre-Curved**
Configurations



11.5F LUMEN
Small lumen design.

TRIPLE LUMEN
Provides reliable access for not only dialysis but also infusion and withdrawal needs.



TAPERED TIP
Reduces vessel trauma and provides an easier insertion.

	CATHETER ONLY	
XTP3114CT	11.5F X 12CM TRI-FLOW®, STRAIGHT	10/BOX
XTP3116CT	11.5F X 15CM TRI-FLOW®, STRAIGHT	10/BOX
XTP3118CT	11.5F X 20CM TRI-FLOW®, STRAIGHT	10/BOX
XTP3119CT	11.5F X 24CM TRI-FLOW®, STRAIGHT	10/BOX
XTP3114CTPC	11.5F X 12CM TRI-FLOW®, PRE-CURVED	10/BOX
XTP3116CTPC	11.5F X 15CM TRI-FLOW®, PRE-CURVED	10/BOX
XTP3118CTPC	11.5F X 20CM TRI-FLOW®, PRE-CURVED	10/BOX

	BASIC SET	
XTP3114MT	11.5F X 12CM TRI-FLOW®, STRAIGHT	5/BOX
XTP3116MT	11.5F X 15CM TRI-FLOW®, STRAIGHT	5/BOX
XTP3118MT	11.5F X 20CM TRI-FLOW®, STRAIGHT	5/BOX
XTP3119MT	11.5F X 24CM TRI-FLOW®, STRAIGHT	5/BOX
XTP3114MTPC	11.5F X 12CM TRI-FLOW®, PRE-CURVED	5/BOX
XTP3116MTPC	11.5F X 15CM TRI-FLOW®, PRE-CURVED	5/BOX
XTP3118MTPC	11.5F X 20CM TRI-FLOW®, PRE-CURVED	5/BOX

	FULL TRAY	
MCDLT3-3500	11.5F X 12CM TRI-FLOW®, STRAIGHT	5/BOX
MCDLT3-4000	11.5F X 15CM TRI-FLOW®, STRAIGHT	5/BOX
MCDLT3-4400	11.5F X 20CM TRI-FLOW®, STRAIGHT	5/BOX
MCDLT3-4424	11.5F X 24CM TRI-FLOW®, STRAIGHT	5/BOX
MCDLT3-3500PC	11.5F X 12CM TRI-FLOW®, PRE-CURVED	5/BOX
MCDLT3-4000PC	11.5F X 15CM TRI-FLOW®, PRE-CURVED	5/BOX
MCDLT3-4400PC	11.5F X 20CM TRI-FLOW®, PRE-CURVED	5/BOX

CATHETER ONLY CONTAINS:

- (1) Catheter
- (1) Dilator
- (3) Injection Caps
- (1) Patient I.D. Card

BASIC SET CONTAINS:

- (1) Catheter
- (1) Introducer Needle
- (1) Guidewire
- (1) Dilator
- (1) Adhesive Wound Dressing
- (3) Injection Caps
- (1) Patient I.D. Card

FULL TRAY CONTAINS:

- (1) Catheter
- (1) Introducer Needle
- (1) Guidewire
- (1) Dilator
- (1) Safety Scalpel
- (1) Adhesive Wound Dressing
- (3) Injection Caps
- (8) 4" x 4" Gauze
- (2) Chlorhexidine
- (2) 5cc Syringe
- (3) Pre-Filled Saline Syringes
- (1) 2-0 Silk Suture
- (3) Safety Hypodermic Needles
- (1) Hemostat
- (1) Fenestrated Drape
- (1) 5cc Lidocaine
- (1) Filter Straw
- (1) Surgical Gloves
- (1) Needle Stick Pad
- (1) Patient I.D. Card