

Coaxial

# Micro

# Stick<sup>®</sup>

Introducer Set

Performance  
By Design

# Coaxial Micro-Stick®


## Medcomp® introduces our new

Through advanced engineering and careful design details, we focused on the system as a whole to create performance by design.

### Micro-Stick® Design

Engineers designed the tip around enhancing the performance of the procedure.

- Longer tip for a more gradual dilation of the vessel.
- New, sleeker tip design features a zero-clearance transition between sheath and dilator, allowing a smoother insertion into the vessel.
- Stiff version features a steel cannula inside the dilator, providing extra rigidity for accessing tough vessels.



*New sleek tip outperformed competitor's structural integrity, ensuring reliable performance you can expect.*



Micro-Stick® is a registered trademark of Medcomp®.

New  
STIFF Version  
Now Available.

## new coaxial Micro-Stick® design.

Medcomp® offers several options, so clinicians can easily tailor their procedure to the patients' needs.

### *Coaxial Micro-Stick®*

- A choice of 4 and 5 French sheath sizes.
- Stiffer version available to ensure access into tough vessels.

### *21 Gauge Super Sharp™ Needle*

- Provides less traumatic insertion.
- Echogenic tip for higher visibility.
- Bevel indicator for confidence in bevel position.
- Echogenic and non-echogenic tips available.

### *Variety of .018 x 45cm Guide Wires*

- Nitinol or Steel wires.
- Radiopaque distal tip enhances visualization.
- Curvature of floppy tip allows flexible ease into the vessel.

Catalog #	Nitinol Wire with Radiopaque Tip and Echogenic Needle
MR190421	4F MICRO-STICK INTRODUCER SET, WITH NITINOL WIRE
MR190521	5F MICRO-STICK INTRODUCER SET, WITH NITINOL WIRE
MR190411	4F MICRO-STICK STIFF INTRODUCER SET, WITH NITINOL WIRE
MR190511	5F MICRO-STICK STIFF INTRODUCER SET, WITH NITINOL WIRE

Catalog #	Steel Wire with Radiopaque Tip and Echogenic Needle
MR190422	4F MICRO-STICK INTRODUCER SET, WITH STEEL WIRE
MR190522	5F MICRO-STICK INTRODUCER SET, WITH STEEL WIRE
MR190412	4F MICRO-STICK STIFF INTRODUCER SET, WITH STEEL WIRE
MR190512	5F MICRO-STICK STIFF INTRODUCER SET, WITH STEEL WIRE

Catalog #	Nitinol Wire with Radiopaque Tip and Non-Echogenic Needle
MR193421	4F MICRO-STICK INTRODUCER SET, WITH 4CM NON-ECHO TIP
MR193521	5F MICRO-STICK INTRODUCER SET, WITH 4CM NON-ECHO TIP
MR193411	4F MICRO-STICK STIFF INTRODUCER SET, WITH 4CM NON-ECHO TIP
MR193511	5F MICRO-STICK STIFF INTRODUCER SET, WITH 4CM NON-ECHO TIP
MR192421	4F MICRO-STICK INTRODUCER SET, WITH 7CM NON-ECHO TIP
MR192521	5F MICRO-STICK INTRODUCER SET, WITH 7CM NON-ECHO TIP
MR192411	4F MICRO-STICK STIFF INTRODUCER SET, WITH 7CM NON-ECHO TIP
MR192511	5F MICRO-STICK STIFF INTRODUCER SET, WITH 7CM NON-ECHO TIP

Catalog #	Steel Wire with Radiopaque tip and Non-Echogenic Needle
MR193422	4F MICRO-STICK INTRODUCER SET, WITH 4CM NON-ECHO TIP
MR193522	5F MICRO-STICK INTRODUCER SET, WITH 4CM NON-ECHO TIP
MR193412	4F MICRO-STICK STIFF INTRODUCER SET, WITH 4CM NON-ECHO TIP
MR193512	5F MICRO-STICK STIFF INTRODUCER SET, WITH 4CM NON-ECHO TIP
MR192422	4F MICRO-STICK INTRODUCER SET, WITH 7CM NON-ECHO TIP
MR192522	5F MICRO-STICK INTRODUCER SET, WITH 7CM NON-ECHO TIP
MR192412	4F MICRO-STICK STIFF INTRODUCER SET, WITH 7CM NON-ECHO TIP
MR192512	5F MICRO-STICK STIFF INTRODUCER SET, WITH 7CM NON-ECHO TIP



**medCOMP®**

1499 Delp Drive  
Harleysville, PA 19438  
Tel: 215-256-4201  
Fax: 215-256-1787  
www.medcompnet.com

CE<sub>0434</sub>  
**ISO13485:2003**  
Certified Quality System