




M
U
L
T
I
C
A
T
H

MULTI CATH

*The Highest Quality of
Innovative Infusion Therapy*



Central Vein
Infusion
Catheters

Multi-Cath

The Highest Quality of Innovative Infusion Therapy

Medcomp® provides product innovation and superior technology to extend the use of catheterization for the treatment and diagnosis of critically ill patients. The Multi-Cath infusion catheter allows the administration of fluids, drugs, and blood products. Double, triple, and quad lumens allow incompatible drugs to be administered through the same catheter. These products are also used for patient monitoring and diagnosis. This supreme catheter is a vital tool for the intensivist where central access is limited but vital.



• *Printed depth markings assists in accurate central venous placement.*

• *Adjustable Wing: Allows anchorage at puncture site regardless of depth of catheter. This minimizes trauma and irritation to the exit site.*



MULTI-CATH

• Medcomp CVC's are available in a variety of lengths, French sizes, and diameters.

• Flexible polyurethane designed to enhance indwelling characteristics.

• Vessel dilators aid in placing the catheter percutaneously thereby ensuring a soft catheter/vein interface.

• New 8.5FR available providing higher flow rates and CT Rated.
(Available 2009)

• Radiopacity allows fluoroscopic confirmation of catheter placement.

MULTI-CATH ORDERING INFORMATION

Single Lumen Catheter Sets (20/Box)

14ga x 15cm Single Lumen	CSP-1406
14ga x 20cm Single Lumen	CSP-1408
14ga x 30cm Single Lumen	CSP-1412
16ga x 15cm Single Lumen	CSP-1606
16ga x 20cm Single Lumen	CSP-1608
16ga x 30cm Single Lumen	CSP-1612
18ga x 15cm Single Lumen	CSP-1806
18ga x 20cm Single Lumen	CSP-1808
18ga x 30cm Single Lumen	CSP-1812
20ga x 13cm Single Lumen	CSP-2005
24ga x 13cm Single Lumen	CSP-2405

Double Lumen Catheter Sets (10/Box)

4FR x 5cm Double Lumen	MC2LP-35S
4FR x 8cm Double Lumen	MC2LP-38S
4FR x 13cm Double Lumen	MC2LP-313S
4FR x 30cm Double Lumen	MC2LP-330S
5FR x 8cm Double Lumen	MC2LP-58S
5FR x 13cm Double Lumen	MC2LP-5S
5FR x 30cm Double Lumen	MC2LP-530S
7FR x 15cm Double Lumen	MC2L-6S
7FR x 20cm Double Lumen	MC2L-8S
7FR x 30cm Double Lumen	MC2L-12S
7.5FR x 15cm Double Lumen	MC2L-715S
7.5FR x 20cm Double Lumen	MC2L-720S
7.5FR x 30cm Double Lumen	MC2L-730S
8FR x 11cm Double Lumen	MC2L8-4S
8FR x 15cm Double Lumen	MC2L8-6S
8FR x 20cm Double Lumen	MC2L8-8S
8FR x 30cm Double Lumen	MC2L8-12S
8.5FR x 11cm Double Lumen	MC2L8.5-4S
8.5FR x 15cm Double Lumen	MC2L8.5-6S
8.5FR x 20cm Double Lumen	MC2L8.5-8S
8.5FR x 30cm Double Lumen	MC2L8.5-12S

Triple Lumen Catheter Sets (10/Box)

5.5FR x 8cm Triple Lumen	MC3LP-58S
5.5FR x 13cm Triple Lumen	MC3LP-5S
5.5FR x 30cm Triple Lumen	MC3LP-530S
7FR x 15cm Triple Lumen	MC3L-6S
7FR x 20cm Triple Lumen	MC3L-8S
7FR x 30cm Triple Lumen	MC3L-12S
7.5FR x 15cm Triple Lumen	MC3L-715S
7.5FR x 20cm Triple Lumen	MC3L-720S
7.5FR x 30cm Triple Lumen	MC3L-730S

Quad Lumen Catheter Sets (10/Box)

8.5 FR x 15cm Quad Lumen	MC4L-6S
8.5 FR x 20cm Quad Lumen	MC4L-8S
8.5 FR x 30cm Quad Lumen	MC4L-12S

MULTI LUMEN SETS CONTAIN:

- (1) Catheter with Clamps
- (1) J/Flex Guidewire
- (1) Introducer Needle
- (2) Injection Caps
 - (3) for Triple Lumen
 - (4) for Quad Lumen
- (1) 5cc Syringe
- (1) Vessel Dilator
- (1) Movable Wing (Double and Triple Lumen Only)

SINGLE LUMEN SETS CONTAIN:

- (1) Catheter with Clamps
- (1) J/Flex Guidewire
- (1) Introducer Needle
- (1) Vessel Dilator
- (1) Movable Wing (14ga - 18ga Only)

Multi Lumen Double Wrapped Full Trays

4FR x 5cm Double	MC2LP-35S-PA
4FR x 8cm Double	MC2LP-38S-PA
5FR x 8cm Double	MC2LP-58S-PA
5FR x 13cm Double	MC2LP-5S-PA
7FR x 20cm Double	MC2L-8S-PA
5FR x 8cm Triple	MC3LP-58S-PA
5.5FR x 13cm Triple	MC3LP-5S-PA
7FR x 20cm Triple	MC3L-8S-PA

MULTI LUMEN DOUBLE WRAPPED FULL STERILE TRAYS CONTAIN:

- (1) Catheter with Clamps
- (1) J/Flex Guidewire (0.035" for Adult, 0.018" for Ped/Neonate)
- (1) Introducer Needle (18G/Adult, 21Gx1-1/2"/Ped, 21x1"/Neonate)
- (2) Injection Caps
 - (3) for Triple Lumen
- (1) Syringe (5cc/Adult, 3cc/Ped, 1cc/Neonate)
- (1) Vessel Dilator
- (1) Suture (3.0/Adult, 4.0 Ped/Neonate) w/needle
- (4) 4" x 4" Gauze
- (1) Needle Holder
- (1) Chloraprep
- (1) 5cc Ampule Lidocaine HCl 1%
- (1) Safety Scalpel
- (1) 22ga x 1-1/2" Needle
- (1) 25ga x 1" Needle
- (1) Fenestrated Drape

Medcomp® offers a full range of Neonatal catheters including micro puncture access.