

# PATIENT SIZE // CATHETER SIZE

## BFS-6 // SMALL PATIENT-RIGHT SIDE


Cuffs are 18CM from distal tip (Arterial),  
21CM from the distal tip (Venous)

## BFR-6 // RIGHT SIDE

Cuffs are 22CM from distal tip (Arterial),  
25CM from the distal tip (Venous)

## BFL-6 // LEFT SIDE OR OBESE PATIENT

Cuffs are 27CM from distal tip (Arterial),  
30CM from the distal tip (Venous)



## BF11V // (BLUE) ADAPTOR (VENOUS)

## BF11A // (RED) ADAPTOR (ARTERIAL)

