

VASCU-PICC®
VASCU-PICC® TL & TLS



VASCU-PICC®
•••••

www.medoompnet.com

INTRODUCING VASCU-PICC® & **NEW** VASCU-PICC® TL & TLS

Ranging in size from 1.9F to 5F, and available with tapered or non-tapered lumens, Medcomp's Vascu-PICC® line of catheters provides a full range of options designed to meet the unique access challenges of patients young and old. Consistent across all designs is a radiopaque polyurethane material, which provides device visibility, high flows and superior catheter strength, ensuring device reliability and patient comfort.



Reduce bleeding at the venotomy site with our tapered design.



Featuring a taperless lumen, ideal for smaller vasculature.



Our smallest, taperless PICC, specifically designed
& FDA indicated for pediatric patients.



FEATURES & BENEFITS

Vascu-PICC®, Vascu-PICC® TL, & Vascu-PICC® TLS

Thermosensitive Material

Softens at body temperature.

Radiopaque Lumen

Identifiable under X-ray.

Easy-to-Read Depth Marks

Easily & accurately trim to the desired length.

Crack Resistant Luers

Ensures patient safety & device longevity.

Polyurethane Material

Provides strength, versatility, & high flow rates.

Non-Tapered Design

Reduces the incidence of vessel occlusion in smaller vasculature.

Note: Non-Tapered design only applicable to the Vascu-PICC® TL & Vascu-PICC® TLS.

Tapered Design

Reduces bleeding at the venotomy site.

Note: Tapered design only applicable to the Vascu-PICC®.



[Ordering Information](#) >



Catalog #	Catheter Only	UOM
MR17013100	3F X 60CM VASCU-PICC®, SINGLE LUMEN	10/BOX
MR17014100	4F X 60CM VASCU-PICC®, SINGLE LUMEN	10/BOX
MR17014200	4F X 60CM VASCU-PICC®, DOUBLE LUMEN	10/BOX
MR17015200	5F X 60CM VASCU-PICC®, DOUBLE LUMEN	10/BOX

Set Contains: (1) Vascu-PICC® w/Stylet (1) Patient Information Pack (1|2) Needle Free Valve(s) (1) Catheter Securement Device (1) Patient Chart Sticker

Catalog #	Basic Set	UOM
MR17013101	3F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX
MR17014101	4F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX
MR17015101	5F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX
MR17014201	4F X 60CM VASCU-PICC®, DOUBLE LUMEN	5/BOX
MR17015201	5F X 60CM VASCU-PICC®, DOUBLE LUMEN	5/BOX
MR17015301	5F X 60CM VASCU-PICC®, TRIPLE LUMEN	5/BOX

Set Contains: (1) Vascu-PICC® w/Stylet (1) Peelable Sheath Introducer (1) Mini Scalpel (1) Guidewire (1) 10cc Syringe (1) Tape Measure (1) Introducer Needle (1|2|3) Needle Free Valve(s) (1) Catheter Securement Device (1) Patient Information Pack (1) Patient Chart Sticker

Catalog #	Long Wire Set	UOM
MR17014102	4F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX
MR17015102	5F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX
MR17015202	5F X 60CM VASCU-PICC®, DOUBLE LUMEN	5/BOX
MR17015302	5F X 60CM VASCU-PICC®, TRIPLE LUMEN	5/BOX

Set Contains: (1) Vascu-PICC® (1) Peelable Sheath Introducer (1) Mini Scalpel (1) Guidewire (1) 10cc Syringe (1) Tape Measure (1) Introducer Needle (1|2|3) Needle Free Valve(s) (1) Catheter Securement Device (1) Patient Information Pack (1) Patient Chart Sticker

Catalog #	Basic Nursing	UOM
MR17013105	3F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX
MR17014105	4F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX
MR17015105	5F X 60CM VASCU-PICC®, SINGLE LUMEN	5/BOX
MR17014205	4F X 60CM VASCU-PICC®, DOUBLE LUMEN	5/BOX
MR17015205	5F X 60CM VASCU-PICC®, DOUBLE LUMEN	5/BOX
MR17015305	5F X 60CM VASCU-PICC®, TRIPLE LUMEN	5/BOX

Set Contains: (1) Vascu-PICC® w/Stylet (1) Peelable Sheath Introducer (1) Safety Scalpel (1) Guidewire (1) 10cc Syringe (1) Tape Measure (1) Safety Introducer Needle (1|2|3) Needle Free Valve(s) (1) Catheter Securement Device (1) Latex-Free Tourniquet (1) Patient Information Pack (1) Patient Chart Sticker

Catalog #	Full Nursing Kit	UOM
MR17013106	3F X 60CM VASCU-PICC®, SINGLE LUMEN	4/BOX
MR17014106	4F X 60CM VASCU-PICC®, SINGLE LUMEN	4/BOX
MR17015106	5F X 60CM VASCU-PICC®, SINGLE LUMEN	4/BOX
MR17014206	4F X 60CM VASCU-PICC®, DOUBLE LUMEN	4/BOX
MR17015206	5F X 60CM VASCU-PICC®, DOUBLE LUMEN	4/BOX
MR17015306	5F X 60CM VASCU-PICC®, TRIPLE LUMEN	4/BOX

Set Contains: (1) Vascu-PICC® w/Stylet (1) Peelable Sheath Introducer (1) Safety Scalpel (1) Guidewire (1) 5cc Syringe (2) Tape Measures (1) Surgical Tape (1) Safety Introducer Needle (1|2|3) Needle Free Valve(s) (1) Catheter Securement Device (1) Latex-Free Tourniquet (3) Pre-Filled Saline Syringes (1) Chlorhexidine (1) Tinted Chlorhexidine (1) 5cc Lidocaine (1) Adhesive Wound Dressing (1) Needle Stick Pad (10) 4" x 4" Gauze (3) Surgical Gloves (1) Filter Straw (1) Scissors (1) Steri-Strips (2) Face Masks (8) 2" x 2" Gauze (1) I.V. Catheter (1) Towel (1) Fenestrated Drape (1) Drape (1) Hypodermic Safety Needle (1) Lidocaine/Saline Stickers (1) Forceps (1) Introducer Needle (1) Patient Information Pack (1) Patient Chart Sticker

Catalog #	Max Barrier Kit	UOM
MR17014106MB	4F X 60CM VASU-PICC®, SINGLE LUMEN	3/BOX
MR17015206MB	5F X 60CM VASCU-PICC®, DOUBLE LUMEN	3/BOX
MR17015306MB	5F X 60CM VASCU-PICC®, TRIPLE LUMEN	3/BOX

Set Contains: (1) Vascu-PICC® w/Stylet (1) Peelable Sheath Introducer (1) Safety Scalpel (1) Safety Introducer Needle (1|2|3) Needle Free Valve(s) (1) Catheter Securement Device (1) Latex-Free Tourniquet (3) Pre-Filled Saline Syringes (1) Tinted Chlorhexidine (1) Chlorhexidine (1) 5cc Lidocaine (1) Adhesive Wound Dressing (1) Needle Stick Pad (10) 4" x 4" Gauze (3) Surgical Gloves (1) Filter Straw (1) Scissors (1) Steri-Strips (2) Face Masks (8) 2" x 2" Gauze (1) I.V. Catheter (1) Towel (1) Fenestrated Drape (1) Drape (1) Hypodermic Safety Needle (1) Lidocaine/Saline Stickers (1) Full Body Drape (1) XL Gown (1) Surgical Face Mask w/ Eye Shield (1) Bouffant Cap (1) Forceps (1) Saline Wipes (1) Introducer Needle (1) Patient Information Pack (1) Patient Chart Sticker (1) 5cc Syringe (1) Probe Cover w/Gel (1) Guidewire (1) Surgical Tape (2) Tape Measures



Catalog #	Catheter Only	UOM
MR81013100	3F X 60CM VASCU-PICC® TL, SINGLE LUMEN	10/BOX
MR81014100	4F X 60CM VASCU-PICC® TL, SINGLE LUMEN	10/BOX
MR81014200	4F X 60CM VASCU-PICC® TL, DOUBLE LUMEN	10/BOX

Set Contains: (1) Vascu-PICC® TL w/Stylet (1) Patient Information Pack (1|2) Needle Free Valve(s) (1) Catheter Securement Device (1) Patient Chart Sticker

Catalog #	Basic Set	UOM
MR81013101	3F X 60CM VASCU-PICC® TL, SINGLE LUMEN	5/BOX
MR81014101	4F X 60CM VASCU-PICC® TL, SINGLE LUMEN	5/BOX
MR81015101	5F X 60CM VASCU-PICC® TL, SINGLE LUMEN	5/BOX
MR81014201	4F X 60CM VASCU-PICC® TL, DOUBLE LUMEN	5/BOX
MR81015201	5F X 60CM VASCU-PICC® TL, DOUBLE LUMEN	5/BOX
MR81015301	5F X 60CM VASCU-PICC® TL, TRIPLE LUMEN	5/BOX

Set Contains: (1) Vascu-PICC® TL w/Stylet (1) Introducer Needle (1) Patient Chart Sticker (1) Peelable Sheath Introducer (1) 10cc Syringe (1) Patient Information Pack (1) Guidewire (1) Mini Scalpel (1) Catheter Securement Device (1|2|3) Needle Free Valve(s) (1) Tape Measure

Catalog #	Basic Nursing	UOM
MR81013105	3F X 60CM VASCU-PICC® TL, SINGLE LUMEN	5/BOX
MR81014105	4F X 60CM VASCU-PICC® TL, SINGLE LUMEN	5/BOX
MR81015105	5F X 60CM VASCU-PICC® TL, SINGLE LUMEN	5/BOX
MR81014205	4F X 60CM VASCU-PICC® TL, DOUBLE LUMEN	5/BOX
MR81015205	5F X 60CM VASCU-PICC® TL, DOUBLE LUMEN	5/BOX
MR81015305	5F X 60CM VASCU-PICC® TL, TRIPLE LUMEN	5/BOX

Set Contains: (1) Vascu-PICC® TL w/Stylet (1) Safety Introducer Needle (1) Tape Measure (1) Peelable Sheath Introducer (1) 10cc Syringe (1) Patient Chart Sticker (1) Guidewire (1) Safety Scalpel (1) Patient Information Guide (1|2|3) Needle Free Valve(s) (1) Latex Free Tourniquet (1) Catheter Securement Device

Catalog #	Max Barrier	UOM
MR81013106MB	3F X 60CM VASCU-PICC® TL, SINGLE LUMEN	3/BOX
MR81014106MB	4F X 60CM VASCU-PICC® TL, SINGLE LUMEN	3/BOX
MR81015106MB	5F X 60CM VASCU-PICC® TL, SINGLE LUMEN	3/BOX
MR81014206MB	4F X 60CM VASCU-PICC® TL, DOUBLE LUMEN	3/BOX
MR81015206MB	5F X 60CM VASCU-PICC® TL, DOUBLE LUMEN	3/BOX
MR81015306MB	5F X 60CM VASCU-PICC® TL, TRIPLE LUMEN	3/BOX

Set Contains: (1) Vascu-PICC® TL w/Stylet (1) 5cc Syringe (1) Adhesive Dressing (1) Peelable Sheath Introducer (1) Safety Scalpel (1) Surgical Tape (1) Guidewire (2) Tape Measure (1) Saline Wipe (3) Pre-Filled Saline Syringe (1) Probe Cover w/Gel (2) Face Masks (1) 5cc Lidocaine (1) Safety Needle (1) Steri-Strips (3) Surgical Gloves (1) Safety Introducer Needle (1) Patient Chart Sticker (1) XL Gown (1) Filter Straw (1) Patient Information Guide (1) Full Body Drape (1) Scissors (1) Lidocaine/Saline Stickers (1) Fenestrated Drape (1) Forceps (1) Needle Stick Pad (1) Drape (1) Surgical Face Mask with Eyeshield (1) Catheter Securement Device (2) Chlorhexidine (1) Bouffant Cap (10) 4" x 4" Gauze (1|2|3) Needle Free Valve(s) (1) Towel (8) 2" x 2" Gauze (1) Introducer Needle (1) Hypodermic Safety Needle (1) Latex Free Tourniquet



Catalog #	Catheter Only	UOM
MR17011100	1.9F X 50CM VASCU-PICC® TLS, SINGLE LUMEN	10/BOX
MR17012600	2.6F X 50CM VASCU-PICC® TLS, DOUBLE LUMEN	10/BOX

Set Contains: (1) Vascu-PICC® TLS w/preloaded Stylet (1) Patient Information Booklet (1|2) Needle Free Valve(s) (1) Patient Chart Sticker (1) Catheter Securement Device

Catalog #	Basic Set	UOM
MR17011121	1.9F X 20CM VASCU-PICC® TLS, SINGLE LUMEN	5/BOX
MR17011101	1.9F X 50CM VASCU-PICC® TLS, SINGLE LUMEN	5/BOX
MR17012621	2.6F X 20CM VASCU-PICC® TLS, DOUBLE LUMEN	5/BOX
MR17012601	2.6F X 50CM VASCU-PICC® TLS, DOUBLE LUMEN	5/BOX

Set Contains: (1) Vascu-PICC® TLS w/preloaded Stylet (1) Tearaway Introducer (1) Latex Free Tourniquet (1) 24" Tape Measure (10) 2" x 2" Gauze (1) Catheter Securement Device (1) 10cc Syringe (1) Patient Chart Sticker (1|2) Needle Free Valve(s)

Catalog #	Procedure Tray	UOM
MR17019000	PROCEDURE TRAY	3/BOX

Set Contains: (3) Pre-Filled Saline Syringes (1) #11 Safety Scalpel (1) Scissors (2) Chlorhexidine (1) Plastic Forceps (1) Surgical Tape (1) 25" x 14" Absorbant Towel (1) 37" x 24" Fenestrated Drape (1) Steri-Strips (1) 36" x 36" CSR Wrap (1) Adhesive Dressing (1) Needle Stick Pad (1) 5cc Lidocaine HCL 1% (1) Filter Straw (1) Lidocaine/Saline Stickers (3) Surgical Gloves (1) Saline Wipes (2) 24" Tape Measures (1) 40" x 70" Drape (2) Face Masks (1) 36" x 24" Drape (1) 5cc Disposable Luer Lock Syringe (1) Bouffant Cap (10) 4" x 4" Gauze (1) Hypodermic Safety Needle (1) 24" x 24" CSR Wrap (8) 2" x 2" Gauze (1) XL Gown

Catalog #	Pediatric Procedure Tray	UOM
MR17019001	PROCEDURE TRAY	3/BOX

Set Contains: (3) Pre-Filled Saline Syringes (1) XL Gown (1) Scissors (2) Chlorhexidine (1) #11 Safety Scalpel (1) Surgical Tape (1) 25" x 14" Absorbant Towel (1) Steel Forceps (1) Steri-Strips (1) 36" x 36" CSR Wrap (1) 19" x 12" Fenestrated Drape (1) Needle Stick Pad (1) 5cc Lidocaine HCL 1% (1) Adhesive Dressing (1) Lidocaine/Saline Stickers (3) Surgical Gloves (1) Filter Straw (2) 24" Tape Measures (1) 5cc Disposable Luer Lock Syringe (1) Saline Wipe (1) 36" x 24" Drape (1) Hypodermic Safety Needle (2) Face Masks (10) 4" x 4" Gauze (1) 24" x 24" CSR Wrap (1) Bouffant Cap (8) 2" x 2" Gauze



1499 Delp Dr.
 Harleysville, PA 19438
 P: 215.256.4201 F: 215.256.1787
 www.medcompnet.com

ISO13485:2003
 Certified Quality System