

CENTRAL VENOUS CATHETER (CVC)
**CHRONOFLEX® 5F
 CATHETER**



Kit Contains

- (1) CT Implantable Port
- (1) Catheter
- (2) Catheter Locks
- (1) Scalpel
- (1) 18ga Introducer Needle
- (1) Vein Pick
- (1) .038" x 70cm "J" Marked Guidewire with Advancer
- (1) Tearaway Introducer
- (1) 10cc Syringe
- (1) Tunneler
- (2) 22ga Huber Needles (1 Straight, 1 Right Angle)
- (1) Blunt Tip Needle
- (1) Port IFU
- (1) Patient Information Pack
- (1) Patient Chart Sticker

Micro-Stick®

- (1) 5F Coaxial Micro-Stick®
- (1) .018" x 45cm Guidewire with Radiopaque Tip
- (1) 21ga Introducer Needle with Echogenic Tip

Small, but Mighty.



	Dignity® CT Port 5F ChronoFlex* Catheter	Dignity® Low Profile CT Port 5F ChronoFlex* Catheter
Attachable Catheter	MRCTI50001	MRCTI50041
Pre-Attached Catheter	MRCTI50801	MRCTI50841

WITH MICRO-STICK® KIT

Attachable Catheter	MICTI50001	MICTI50041
Pre-Attached Catheter	MICTI50801	MICTI50841

WITH SILICONE FILLED SUTURE HOLES

Attachable Catheter	MRCTI5000S	MRCTI5004S
Pre-Attached Catheter	MRCTI5080S	MRCTI5084S

WITH SILICONE FILLED SUTURE HOLES AND MICRO-STICK® KIT

Attachable Catheter	MICTI5000S	MICTI5004S
Pre-Attached Catheter	MICTI5080S	MICTI5084S



	Dignity® Mini CT Port 5F ChronoFlex* Catheter	Dignity® Mini CT Port 5F ChronoFlex* Catheter
Attachable Catheter	MRCTI50041M	MRCTI5004SM
Pre-Attached Catheter	MRCTI50841M	MRCTI5084SM

WITH MICRO-STICK® KIT AND SILICONE FILLED SUTURE HOLES

Attachable Catheter	MICTI50041M	MICTI5004SM
Pre-Attached Catheter	MICTI50841M	MICTI5084SM



	ProFuse® CT Port 5F ChronoFlex* Catheter	ProFuse® CT Port 5F ChronoFlex* Catheter
Attachable Catheter	MRCTT50041	MRCTT5004S
Pre-Attached Catheter	MRCTT50841	MRCTT5084S

WITH MICRO-STICK® KIT AND SILICONE FILLED SUTURE HOLES

Attachable Catheter	MICTT50041	MICTT5004S
Pre-Attached Catheter	MICTT50841	MICTT5084S

* CHRONOFLEX® AL85AB20 (20% Barium) – main ingredient in the catheter lumen.
 CHRONOFLEX® is a registered trademark of AdvanSource Biomaterials, Inc.

CENTRAL VENOUS CATHETER (CVC)
**CHRONOFLEX® 5F
 CATHETER**



CHRONOFLEX® 5F CATHETER

You need enhanced performance and low complication rates. Your patients need comfort.

The ChronoFlex® 5F catheter is the solution. The smallest, most versatile CT power injectable catheter on the market, the 5F can be used with our entire family of Medcomp® CT ports. Ordering is simple with our full range of vascular access products and implantable ports.

Advanced


When it comes to vascular access, smaller is almost always better. The ChronoFlex® 5F catheter offers reliable functionality that both patients and clinicians can depend on.

Versatile

This catheter can be used for chest or arm placement and is available with Medcomp's family of CT implantable ports.

Uncomplicated

The smallest of the power injectable kind, the 5F takes up less space in the vein which improves venous flow and reduces the risk of complications.



Our smallest innovation yet.

