

# TRI-FLOW™ CATHETERS

## FEATURES & BENEFITS

Bio-Flex™ material reduces kinking due to its excellent memory.

Integral soft tip reduces occurrence of stenosis and erosion.

Suture wing allows secure fixation.

Third lumen gives greater flexibility.

Priming volume designations for clear identification.

Length and size designations for clear identification.

## Tri-Flow™

Now enjoy triple lumen versatility or double lumen reliability with our Tri-Flow™ Catheter. The two lumens are ideal for all continuous therapies. The third lumen gives greater flexibility for fluid administration, vasoactive drug therapy, nutrition and blood sampling. The softness minimizes kinking in femoral punctures and limits stenosis/thrombus formation in subclavian applications.

## ORDERING INFORMATION

<i>Tri-Flow™</i>	CATHETERS (10/BOX)	SETS (5/BOX)	TRAYS (5/BOX)
11.5FR x 12cm (4 <sup>3/4</sup> " Triple Lumen Catheter	XTP3114CT	XTP3114MT	MCDLT3-3500
11.5FR x 15cm (6" Triple Lumen Catheter	XTP3116CT	XTP3116MT	MCDLT3-4000
11.5FR x 20cm (8" Triple Lumen Catheter	XTP3118CT	XTP3118MT	MCDLT3-4400

*\*Also Available in Customized Trays*

### CATHETERS CONTAIN:

(1) Catheter (3) Injection Caps (1) Vessel Dilator

### SETS CONTAIN:


(1) Catheter (3) Injection Caps (1) Vessel Dilator (1) J/Flex Guidewire (1) Introducer Needle (1) Adhesive Wound Dressing

### TRAYS CONTAIN:

(1) Catheter (3) Injection Caps (1) Vessel Dilator (1) J/Flex Guidewire (1) Introducer Needle (1) Scalpel (1) Adhesive Wound Dressing (1) Silk Suture (2) 5cc Syringes (2) 10cc Syringes (1) Needle Holder (3) Hypodermic Needles (1) Ampule Lidocaine (1) Pair Surgical gloves (4) Gauze Sponges (3) Povidone Iodine Swabsticks (1) Povidone Iodine Ointment (1) Fenestrated Drape (2) Sterile Wraps



## APPROVED FOR DIALYSIS



Tri-Flow™ Catheter 11.5FR x 15cm (6")

TRI-FLOW™