



FREE FLOW™



No Sideholes / Open Tips



FREE-FLOW™

Short & Long Term Hemodialysis Catheter

*As innovation continues
superior products come to life.*

Medcomp® presents its premier *open tip* dialysis catheter, the Free Flow™. Through the use of leading technology and unprecedented engineering the industry's first true 360° coaxial flow design has come to life, eliminating the need for side holes and reversing of catheter lines due to arterial insufficiency. The *lack of sideholes* preserves catheter locking solution concentration and minimizes locking solution leaching. This unique design increases the overall patency and performance of the catheter.



Raulerson IJ

Soft Tip Stylet (Short Term)

Integral soft tip reduces the risk of vessel damage.



FREE FLOW™ ORDERING INFORMATION

Long-Term Free Flow

15.5FR Straight	Sets (5/Box)
15.5FR x 24cm Straight Long-Term Free Flow Catheter with Stylet Set	MFFL1524S
15.5FR x 28cm Straight Long-Term Free Flow Catheter with Stylet Set	MFFL1528S
15.5FR x 32cm Straight Long-Term Free Flow Catheter with Stylet Set	MFFL1532S
15.5FR x 36cm Straight Long-Term Free Flow Catheter with Stylet Set	MFFL1536S

15.5FR Pre-Curved	Sets (5/Box)
15.5FR x 24cm Pre-Curved Long-Term Free Flow Catheter with Stylet Set	MFFL1524PC
15.5FR x 28cm Pre-Curved Long-Term Free Flow Catheter with Stylet Set	MFFL1528PC
15.5FR x 32cm Pre-Curved Long-Term Free Flow Catheter with Stylet Set	MFFL1532PC
15.5FR x 36cm Pre-Curved Long-Term Free Flow Catheter with Stylet Set	MFFL1536PC

Set contains: (1) Catheter with Stylet, (1) Grasper Tuner, (1) Peel-Away Sheath, (2) Vessel Dilators, (1) J-Flex Guidewire, (1) Introducer Needle, (2) Injection Caps

Short-Term Free Flow

12.5FR Straight	Sets (5/Box)
12.5FR X 15cm Straight Short-Term Free Flow with Stylet Catheter Set	MFSS1215S
12.5FR X 20cm Straight Short-Term Free Flow with Stylet Catheter Set	MFSS1220S
12.5FR X 24cm Straight Short-Term Free Flow with Stylet Catheter Set	MFSS1224S

15.5FR Straight	Sets (5/Box)
15.5FR X 15cm Straight Short-Term Free Flow with Stylet Catheter Set	MFSS1515S
15.5FR X 20cm Straight Short-Term Free Flow with Stylet Catheter Set	MFSS1520S
15.5FR X 24cm Straight Short-Term Free Flow with Stylet Catheter Set	MFSS1524S

12.5FR Raulerson IJ	Sets (5/Box)
12.5FR X 12cm IJ Short-Term Free Flow with Stylet Catheter Set	MFSS1212IJ
12.5FR X 15cm IJ Short-Term Free Flow with Stylet Catheter Set	MFSS1215IJ
12.5FR X 20cm IJ Short-Term Free Flow with Stylet Catheter Set	MFSS1220IJ

15.5FR Raulerson IJ	Sets (5/Box)
15.5FR X 12cm IJ Short-Term Free Flow with Stylet Catheter Set	MFSS1512IJ
15.5FR X 15cm IJ Short-Term Free Flow with Stylet Catheter Set	MFSS1515IJ
15.5FR X 20cm IJ Short-Term Free Flow with Stylet Catheter Set	MFSS1520IJ

12.5FR IJ with 2 Wings	Sets (5/Box)
12.5FR x 12cm IJ Short-Term Free Flow with Stylet and 2 Suture Wings	MFSS1212IJ-2
12.5FR x 15cm IJ Short-Term Free Flow with Stylet and 2 Suture Wings	MFSS1215IJ-2
12.5FR x 20cm IJ Short-Term Free Flow with Stylet and 2 Suture Wings	MFSS1220IJ-2

15.5FR IJ with 2 Wings	Sets (5/Box)
15.5FR x 12cm IJ Short-Term Free Flow with Stylet and 2 Suture Wings	MFSS1512IJ-2
15.5FR x 15cm IJ Short-Term Free Flow with Stylet and 2 Suture Wings	MFSS1515IJ-2
15.5FR x 20cm IJ Short-Term Free Flow with Stylet and 2 Suture Wings	MFSS1520IJ-2

Set contains: (1) Catheter with Stylet, (3) Vessel Dilators, (1) J-Flex Guidewire, (1) Introducer Needle, (2) Injection Caps