



Introducing:

The SLX Catheter



SILICONE MATERIAL

SAME TRUSTED

14^{FR}
DIAMETER

AVAILABLE IN

30^{CM}
LENGTH

CAPABLE OF

500^{ML}
MIN
FLOW RATES

Direct from Production to Procedure.

www.medcompnet.com



Ordering Information:

Catalog #	Catheter Only	
MC061442	14F x 15CM SLX DOUBLE LUMEN, CATHETER ONLY	10/Box
MC061443	14F x 20CM SLX DOUBLE LUMEN, CATHETER ONLY	10/Box
MC061444	14F x 24CM SLX DOUBLE LUMEN, CATHETER ONLY	10/Box
MC061445	14F x 30CM SLX DOUBLE LUMEN, CATHETER ONLY	10/Box

Catalog #	Basic Set	
MC061402	14F x 15CM SLX DOUBLE LUMEN, BASIC SET	5/Box
MC061403	14F x 20CM SLX DOUBLE LUMEN, BASIC SET	5/Box
MC061404	14F x 24CM SLX DOUBLE LUMEN, BASIC SET	5/Box
MC061405	14F x 30CM SLX DOUBLE LUMEN, BASIC SET	5/Box

Catalog #	Full Tray	
MC061422	14F x 15CM SLX DOUBLE LUMEN, FULL TRAY	5/Box
MC061423	14F x 20CM SLX DOUBLE LUMEN, FULL TRAY	5/Box
MC061424	14F x 24CM SLX DOUBLE LUMEN, FULL TRAY	5/Box
MC061425	14F x 30CM SLX DOUBLE LUMEN, FULL TRAY	5/Box

CATHETER ONLY CONTAINS:

- (1) Catheter w/ Stylet
- (2) Injection Caps (1) Catheter Clip

BASIC SET CONTAINS:

- (1) Catheter w/ Stylet (1) Introducer Needle (1) Guidewire (2) Dilators (1) Scalpel (1) Tearaway Sheath Introducer (1) Adhesive Wound Dressing (2) Injection Caps (1) Catheter Clip (1) 2-0 Suture

FULL TRAY CONTAINS:

- (1) Catheter w/ Stylet (1) Safety Introducer Needle (1) Guidewire (2) Dilators (1) Safety Scalpel (1) Tearaway Sheath Introducer (1) Adhesive Wound Dressing (2) Injection Caps (1) Catheter Clip (4) 4"x4" Gauze (2) Chlorhexidine (1) 5cc Syringe (2) 10cc Pre-Filled Syringes (1) 2-0 Silk Suture (3) Safety Hypodermic Needles (1) Hemostat (1) Fenestrated Drape (1) 5cc Lidocaine (1) Filter Straw (1) Surgical Gloves

FEATURES & BENEFITS:

- Larger 14F luminal design allows for improved flow rates.
- Silicone material provides biocompatibility and kink resistance.
- Fast Track™ Stylet and round catheter permits easier insertion and less patient discomfort.



Tel: 215.256.4201 Fax: 215-256-1787

www.medcompnet.com