



SPLIT CATH® III

Indications for Use: The Medcomp® Split Cath® III is indicated for use in attaining Long-Term vascular access for hemodialysis and apheresis in the adult population. It may be inserted percutaneously and is primarily placed in the internal jugular vein. Alternate insertion sites include the subclavian vein and inferior vena cava as required. Catheters greater than 40cm are intended for femoral vein or inferior vena cava insertion. Translumbar insertion via inferior vena cava is indicated when all other access sites are identified as non-viable.

SPLIT CATH® WITHOUT SIDHOLES

Indications for Use: The Medcomp® Split Cath® III is indicated for use in attaining Long-Term vascular access for Hemodialysis and Apheresis. It may be inserted percutaneously and is primarily placed in the internal jugular vein. Alternate insertion sites include the subclavian vein. Catheters greater than 40cm are intended for femoral vein insertion.

SPLIT CATH® III, SPLIT CATH® WITHOUT SIDHOLES

Contraindications: This catheter is intended for Long-Term vascular access only and should not be used for any purpose other than indicated in these instructions.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.



Polyurethane
Material

Split Tip
Tip Design

10F, 14F, 16F
French Sizes

**Straight,
Pre-Curved**
Configurations

Antegrade
Insertion Type



DOUBLE "D" INTERNAL LUMEN DESIGN



ARTERIOVENOUS STYLET WEAVE



HOLLOW, FAST TRACK STYLET



SPLITTABLE LUMENS



Refer to the Table of Contents for Important Risk Information regarding this device.

Long Term



STRAIGHT		
ASPC20-3E.	14F X 20CM SPLIT CATH® III	5/BOX
ASPC22-3E.	14F X 22CM SPLIT CATH® III	5/BOX
ASPC24-3E.	14F X 24CM SPLIT CATH® III	5/BOX
ASPC28-3E.	14F X 28CM SPLIT CATH® III	5/BOX
ASPC32-3E.	14F X 32CM SPLIT CATH® III	5/BOX
ASPC36-3E.	14F X 36CM SPLIT CATH® III	5/BOX
ASPC40-3E.	14F X 40CM SPLIT CATH® III	5/BOX
ASPC55-3E.	14F X 55CM SPLIT CATH® III	5/BOX
ASPC2416-3E.	16F X 24CM SPLIT CATH® III	5/BOX
ASPC2816-3E.	16F X 28CM SPLIT CATH® III	5/BOX
ASPC3216-3E.	16F X 32CM SPLIT CATH® III	5/BOX
ASPC3616-3E.	16F X 36CM SPLIT CATH® III	5/BOX
ASPC4016-3E.	16F X 40CM SPLIT CATH® III	5/BOX
ASPC5516-3E.	16F X 55CM SPLIT CATH® III	5/BOX

PRE-CURVED		
ASPC24-3PCE.	14F X 24CM PC SPLIT CATH® III	5/BOX
ASPC28-3PCE.	14F X 28CM PC SPLIT CATH® III	5/BOX
ASPC32-3PCE.	14F X 32CM PC SPLIT CATH® III	5/BOX
ASPC36-3PCE.	14F X 36CM PC SPLIT CATH® III	5/BOX
ASPC2416-3PCE.	16F X 24CM PC SPLIT CATH® III	5/BOX
ASPC2816-3PCE.	16F X 28CM PC SPLIT CATH® III	5/BOX
ASPC3216-3PCE.	16F X 32CM PC SPLIT CATH® III	5/BOX
ASPC3616-3PCE.	16F X 36CM PC SPLIT CATH® III	5/BOX

STRAIGHT WITHOUT SIDEHOLES		
ASPC24-3WOE.	14F X 24CM SPLIT CATH® III W/O SIDEHOLES	5/BOX
ASPC28-3WOE.	14F X 28CM SPLIT CATH® III W/O SIDEHOLES	5/BOX
ASPC32-3WOE.	14F X 32CM SPLIT CATH® III W/O SIDEHOLES	5/BOX
ASPC36-3WOE.	14F X 36CM SPLIT CATH® III W/O SIDEHOLES	5/BOX
ASPC40-3WOE.	14F X 40CM SPLIT CATH® III W/O SIDEHOLES	5/BOX
ASPC55-3WOE.	14F X 55CM SPLIT CATH® III W/O SIDEHOLES	5/BOX

PRE-CURVED WITHOUT SIDEHOLES		
ASPC24-3PCWOE.	14F X 24CM PC SPLIT CATH® III W/O SIDEHOLES	5/BOX
ASPC28-3PCWOE.	14F X 28CM PC SPLIT CATH® III W/O SIDEHOLES	5/BOX
ASPC32-3PCWOE.	14F X 32CM PC SPLIT CATH® III W/O SIDEHOLES	5/BOX
ASPC36-3PCWOE.	14F X 36CM PC SPLIT CATH® III W/O SIDEHOLES	5/BOX

STRAIGHT WITH DUAL STYLETS		
ASPC24-3DSE.	14F X 24CM SPLIT CATH® III W/ DUAL STYLET	5/BOX
ASPC28-3DSE.	14F X 28CM SPLIT CATH® III W/ DUAL STYLET	5/BOX
ASPC32-3DSE.	14F X 32CM SPLIT CATH® III W/ DUAL STYLET	5/BOX
ASPC36-3DSE.	14F X 36CM SPLIT CATH® III W/ DUAL STYLET	5/BOX
ASPC40-3DSE.	14F X 40CM SPLIT CATH® III W/ DUAL STYLET	5/BOX
ASPC55-3DSE.	14F X 55CM SPLIT CATH® III W/ DUAL STYLET	5/BOX

STRAIGHT 10F		
ASPC15P-XL	10F X 15CM SPLIT CATH® III	5/BOX
ASPC18P-XL	10F X 18CM SPLIT CATH® III	5/BOX
ASPC24P-XL	10F X 24CM SPLIT CATH® III	5/BOX

TRANSLUMBAR		
ASPC40-3TLE.	14F X 40CM TRANSLUMBAR KIT	5/BOX
ASPC55-3TLE.	14F X 55CM TRANSLUMBAR KIT	5/BOX

TUNNELER		
AST	SPLIT CATH® III TUNNELER	5/BOX
ASTSG	SPLIT CATH® III SURE GRIP TUNNELER	5/BOX

STRAIGHT & PRE-CURVED BASIC SET CONTAINS:

- (1) Catheter
- (1) Introducer Needle
- (1/2) Guidewire(s)
- (2) Dilators
- (1) Tunneler
- (1) Valved Peelable Introducer
- (2) End Caps

STRAIGHT & PRE-CURVED WITHOUT SIDEHOLES

- SET CONTAINS:
- (1) Catheter
 - (1) Introducer Needle
 - (1/2) Guidewire(s)
 - (2) Dilators
 - (1) Tunneler
 - (1) Valved Peelable Introducer
 - (2) End Caps

STRAIGHT WITH DUAL STYLETS SET CONTAINS:

- (1) Catheter w/ Stylets
- (1) Introducer Needle
- (1) Guidewire
- (2) Dilators
- (1) Tunneler
- (1) Valved Peelable Introducer
- (2) End Caps

STRAIGHT 10F BASIC SET CONTAINS:

- (1) Catheter
- (1) Introducer Needle
- (1) Guidewire
- (2) Dilators
- (1) Scalpel
- (1) Tunneler
- (1) Peelable Introducer
- (1) Adhesive Wound Dressing
- (2) End Caps

TRANSLUMBAR KIT CONTAINS:

- (1) Catheter
- (3) Dilators
- (2) End Caps
- (1) Scalpel
- (1) Peelable Introducer
- (1) Introducer Needle
- (2) Guidewires
- (1) Adhesive Wound Dressing
- (1) Tunneler
- (1) Stiffening Stylet