

SYMETREX[®]

LONG TERM HEMODIALYSIS CATHETER

PATIENT INFORMATION CARD

(Please show to Dialysis Provider)

GREEN MEANS GO[™]

SYMETREX[®] Catheters are REVERSIBLE*

<1% Recirculation in Forward and Reverse Flow

**when tested in vitro, in vitro data on file*

Reverse lines as needed

Patient Name: _____ Insertion Date: ____/____/____

SYMETREX® Long Term Hemodialysis Catheter Lot #: _____

Physician: _____ Phone #: _____

Warning: Do not use acetone on any part of the catheter tubing. Chlorhexedine is the preferred alternative. Do not use alcohol to lock catheter because alcohol has been shown to degrade polyurethane catheter tubing with prolonged exposure.

Consult Instructions for Use

 1499 Delp Drive, Harleysville, PA 19438, USA PN2627 Rev. E 11/18

www.medcompnet.com