



LOCKING DILATOR
Assures stability of catheter during insertion.

MULTIPLE CONFIGURATIONS
Available in two lengths and French sizes.

ZERO-CLEARANCE SHEATH/DILATOR TRANSITION
Facilitates a less traumatic insertion.

BASIC SET		
MC5CAVH3	5F X 7.6CM C.A.V.H CATHETER, STRAIGHT	5/BOX
MC8CAVH4F	8F X 11CM C.A.V.H. CATHETER, STRAIGHT	5/BOX
MC8CAVH6	8F X 15CM C.A.V.H. CATHETER, STRAIGHT	5/BOX

FULL TRAY		
MCT-106	8F X 15CM C.A.V.H. CATHETER, STRAIGHT	5/BOX

- BASIC SET CONTAINS:**
- (2) Catheters w/ Dilators
 - (2) Injection Caps
 - (2) Introducer Needles
 - (2) Guidewires
 - (1) Scalpel
 - (2) Adhesive Wound Dressings
 - (1) Patient I.D. Card

- FULL TRAY CONTAINS:**
- (2) Catheters w/ Dilators
 - (2) Injection Caps
 - (2) Introducer Needles
 - (2) Guidewires
 - (1) Scalpel
 - (1) Adhesive Wound Dressing
 - (1) Suture
 - (3) 10cc Syringes
 - (2) Hypodermic Needles
 - (2) 5cc Lidocaine
 - (1) Surgical Gloves
 - (12) 4" x 4" Gauze
 - (1) Povidone Iodine Solution
 - (1) Fenestrated Drape
 - (2) Sterile Wraps
 - (2) Prep Sponge Applicators
 - (4) Towels
 - (1) Patient I.D. Card