



DOUBLE "D" INTERNAL LUMEN DESIGN

Polyurethane  
Material

Split Tip  
Tip Design

13F  
French Size

Straight,  
Pre-Curved,  
Raulerson IJ  
Configurations




**SPLIT TIP DESIGN**  
Improves patency, flow and lowers recirculation.



**ARTERIOVENOUS WEAVE WITH DUAL STYLET**  
Removes the incidence of buckling during insertion.



**RAULERSON IJ CONFIGURATION**  
Reduces infection and improves patient comfort during treatment.



	CATHETER ONLY	
DSP134C	13F X 12CM DUO-SPLIT®, STRAIGHT	10/BOX
DSP136C	13F X 15CM DUO-SPLIT®, STRAIGHT	10/BOX
DSP138C	13F X 20CM DUO-SPLIT®, STRAIGHT	10/BOX
DSP139C	13F X 24CM DUO-SPLIT®, STRAIGHT	10/BOX
DSP130C	13F X 30CM DUO-SPLIT®, STRAIGHT	10/BOX
DSP134PC	13F X 12CM DUO-SPLIT®, PRE-CURVED	10/BOX
DSP136PC	13F X 15CM DUO-SPLIT®, PRE-CURVED	10/BOX
DSP138PC	13F X 20CM DUO-SPLIT®, PRE-CURVED	10/BOX
DSP139PC	13F X 24CM DUO-SPLIT®, PRE-CURVED	10/BOX
DSP134JC	13F X 12CM DUO-SPLIT®, RAULERSON IJ	10/BOX
DSP136JC	13F X 15CM DUO-SPLIT®, RAULERSON IJ	10/BOX
DSP138JC	13F X 20CM DUO-SPLIT®, RAULERSON IJ	10/BOX
DSP139JC	13F X 24CM DUO-SPLIT®, RAULERSON IJ	10/BOX

	FULL TRAY	
DSP134T	13F X 12CM DUO-SPLIT®, STRAIGHT	5/BOX
DSP136T	13F X 15CM DUO-SPLIT®, STRAIGHT	5/BOX
DSP138T	13F X 20CM DUO-SPLIT®, STRAIGHT	5/BOX
DSP139T	13F X 24CM DUO-SPLIT®, STRAIGHT	5/BOX
DSP130T	13F X 30CM DUO-SPLIT®, STRAIGHT	5/BOX
DSP134PT	13F X 12CM DUO-SPLIT®, PRE-CURVED	5/BOX
DSP136PT	13F X 15CM DUO-SPLIT®, PRE-CURVED	5/BOX
DSP138PT	13F X 20CM DUO-SPLIT®, PRE-CURVED	5/BOX
DSP139PT	13F X 24CM DUO-SPLIT®, PRE-CURVED	5/BOX
DSP134JT	13F X 12CM DUO-SPLIT®, RAULERSON IJ	5/BOX
DSP136JT	13F X 15CM DUO-SPLIT®, RAULERSON IJ	5/BOX
DSP138JT	13F X 20CM DUO-SPLIT®, RAULERSON IJ	5/BOX
DSP139JT	13F X 24CM DUO-SPLIT®, RAULERSON IJ	5/BOX

	BASIC SET	
DSP134S	13F X 12CM DUO-SPLIT®, STRAIGHT	5/BOX
DSP136S	13F X 15CM DUO-SPLIT®, STRAIGHT	5/BOX
DSP138S	13F X 20CM DUO-SPLIT®, STRAIGHT	5/BOX
DSP139S	13F X 24CM DUO-SPLIT®, STRAIGHT	5/BOX
DSP130S	13F X 30CM DUO-SPLIT®, STRAIGHT	5/BOX
DSP134PS	13F X 12CM DUO-SPLIT®, PRE-CURVED	5/BOX
DSP136PS	13F X 15CM DUO-SPLIT®, PRE-CURVED	5/BOX
DSP138PS	13F X 20CM DUO-SPLIT®, PRE-CURVED	5/BOX
DSP139PS	13F X 24CM DUO-SPLIT®, PRE-CURVED	5/BOX
DSP134JS	13F X 12CM DUO-SPLIT®, RAULERSON IJ	5/BOX
DSP136JS	13F X 15CM DUO-SPLIT®, RAULERSON IJ	5/BOX
DSP138JS	13F X 20CM DUO-SPLIT®, RAULERSON IJ	5/BOX
DSP139JS	13F X 24CM DUO-SPLIT®, RAULERSON IJ	5/BOX

CATHETER ONLY CONTAINS:  
(1) Catheter w/ Dual Stylet  
(2) Injection Caps  
(2) Dilators  
(1) Patient I.D. Card

BASIC SET CONTAINS:  
(1) Catheter w/ Dual Stylet  
(2) Injection Caps  
(2) Dilators  
(1) Guidewire  
(1) Introducer Needle  
(1) Adhesive Wound Dressing  
(1) Patient I.D. Card

FULL TRAY CONTAINS:  
(1) Catheter w/ Dual Stylet  
(2) Injection Caps  
(2) Dilators  
(1) Guidewire  
(1) Safety Introducer Needle  
(1) Safety Scalpel  
(1) Adhesive Wound Dressing  
(1) 2-0 Silk Suture  
(1) 5cc Syringe  
(1) Pre-Filled Saline Syringe  
(3) Safety Hypodermic Needles  
(1) 5cc Lidocaine  
(1) Filter Straw  
(1) Surgical Gloves  
(8) 4" x 4" Gauze  
(1) Chlorhexidine  
(1) Fenestrated Drape  
(1) Hemostat  
(1) Needle Stick Pad  
(1) Patient I.D. Card