

Midline Ordering Codes

CATALOG #	MIDLINE CATH ONLY - CT VERSION	QTY/BOX
MLCTD4S200	4F x 20CM CT MIDLINE, SINGLE LUMEN	10
MLCTD5D200	5F x 20CM CT MIDLINE, DUAL LUMEN	10

CATH ONLY

- (1) Midline w/Stylet (1) Catheter Securement Device
- (1) Adaptor w/Sideport (1|2) Needle Free Valve(s)

CATALOG #	MIDLINE BASIC SET - NON CT	QTY/BOX
ML3S20	3F x 20CM MIDLINE, SINGLE LUMEN	5
ML4S20	4F x 20CM MIDLINE, SINGLE LUMEN	5
ML5S20	5F x 20CM MIDLINE, SINGLE LUMEN	5
ML4D20	4F x 20CM MIDLINE, DUAL LUMEN	5
ML5D20	5F x 20CM MIDLINE, DUAL LUMEN	5

CATALOG #	MIDLINE BASIC SET - CT VERSION	QTY/BOX
MLCTD4S20	4F x 20CM CT MIDLINE, SINGLE LUMEN	5
MLCTD5D20	5F x 20CM CT MIDLINE, DUAL LUMEN	5

BASIC SET

- (1) Midline w/ Stylet (1) 10cc Syringe (1) Catheter Securement Device
- (1) Adaptor w/ Sideport (1) Guidewire (1) Safety Scalpel
- (1) Peelable Introducer (1|2) Needle Free Valve(s) (1) Tape Measure
- (1) Safety Introducer Needle

CATALOG #	MIDLINE FULL NURSING TRAY - NON CT	QTY/BOX
ML3S20FN	3F x 20CM MIDLINE, SINGLE LUMEN	5
ML4S20FN	4F x 20CM MIDLINE, SINGLE LUMEN	5
ML5S20FN	5F x 20CM MIDLINE, SINGLE LUMEN	5
ML4D20FN	4F x 20CM MIDLINE, DUAL LUMEN	5
ML5D20FN	5F x 20CM MIDLINE, DUAL LUMEN	5

CATALOG #	MIDLINE FULL NURSING TRAY - CT VERSION	QTY/BOX
MLCT4S20FN	4F x 20CM CT MIDLINE, SINGLE LUMEN	5
MLCT5D20FN	5F x 20CM CT MIDLINE, DUAL LUMEN	5

FULL NURSING TRAY

- (1) Midline w/ Stylet (3) 4" x 4" Gauze (1) Safety Scalpel (1) 20" X 30" Drape
- (1) Adaptor w/ Sideport (3) 2" x 2" Gauze (2) Tape Measures (2) Face Masks
- (1) Tourniquet (1) 10cc Syringe (1) Surgical Tape (1) Surgical Gloves
- (1) Peelable Introducer (2) Pre-Filled Saline Syringes (1) Chlorhexidine (1) Scissors
- (1) Introducer Needle (1) Guidewire (1) 5cc Lidocaine (1) Adhesive Wound Dressing
- (1) Safety Introducer Needle (1|2) Needle Free Valve(s) (1) Filter Straw (1) Needle Stick Pad
- (1) Safety Hypodermic Needle (1) Catheter Securement Device (1) Fenestrated Drape (1) SteriStrips

CATALOG #	MIDLINE MAX BARRIER TRAY - NON CT	QTY/BOX
ML4S20MB	4F x 20CM MIDLINE, SINGLE LUMEN	3
ML5S20MB	5F x 20CM MIDLINE, SINGLE LUMEN	3
ML4D20MB	4F x 20CM MIDLINE, DUAL LUMEN	3
ML5D20MB	5F x 20CM MIDLINE, DUAL LUMEN	3

CATALOG #	MIDLINE MAX BARRIER TRAY - CT VERSION	QTY/BOX
MLCT4S20MB	4F x 20CM CT MIDLINE, SINGLE LUMEN	3
MLCT5D20MB	5F x 20CM CT MIDLINE, DUAL LUMEN	3

MAX BARRIER TRAY

- (1) Midline w/ Stylet (1) Guidewire (1) 53" X 77" Max Barrier Drape
- (1) Adaptor w/ Sideport (1|2) Needle Free Valve(s) (1) Face Mask
- (1) Tourniquet (1) Catheter Securement Device (1) Face Mask with EYESHIELD
- (1) Peelable Introducer (1) Safety Scalpel (1) Probe Cover w/ Gel
- (1) Introducer Needle (2) Tape Measures (1) XL Gown
- (1) Safety Introducer Needle (1) Surgical Tape (1) Bouffant Cap
- (1) Safety Hypodermic Needle (2) Chlorhexidine (1) Scissors
- (3) 4" x 4" Gauze (1) 5cc Lidocaine (1) Adhesive Wound Dressing
- (3) 2" x 2" Gauze (1) Fenestrated Drape (1) Needle Stick Pad
- (1) 10cc Syringe (1) 20" X 30" Drape (1) SteriStrips
- (2) Pre-Filled Saline Syringes

MIDLINE


Short -Term Access



MIDLINE


The **Medcomp® Midline** provides short-term (1-4 weeks) peripheral access to the venous system for selected intravenous therapies and blood sampling. The use of **Midlines** provides an access solution with lower infection and complication rates in comparison to PIVs and CVCs⁽¹⁻³⁾. In addition, **Midlines** do not require an X-ray after placement, can be inserted when a patient's vasculature does not allow a PICC to be completely inserted and also have the potential to reduce needle sticks. Available in multiple packaged configurations, including Max Barrier, the **Medcomp® Midline** is a great solution for short-term access.

Short -Term Access



CT Midline with MST Placement

"Midline" labeling on I.D. Tags



Markings for catheter trimming

Thermosensitive Polyurethane material



Printed identification and purple color denote maximum flow rates for contrast injection.

1. N. Richard Anderson, RN., CRNI; Journal of Infusion Nursing (Vol. 25, No. 5, September/October 2004)
 2. Mayo Clinic Proceedings (September 2006; 81 (9): 1159-1171)
 3. Midlines: An Essential Tool in CLABSI Reduction by: Robert B. Dawson, MSA, BSN, RN, CRNI, CPU, VA-BC, and Nancy L. Moureau, BSN, RN, CRNI, CPU, VA-BC, (March 2013) www.infectioncontrolday.com