

DUO-SPLIT[®]

Short-Term Hemodialysis Catheter



FEATURES & BENEFITS

LUMEN DESIGN

Double "D" Internal



TIP DESIGN

Split Tip



AVAILABLE IN

13F



CONFIGURATIONS

Straight



Pre-Curved



Raulerson IJ



DUO-SPLIT®

Short-Term Hemodialysis Catheter

CATHETER ONLY

DSP134C	13F X 12CM DUO SPLIT® CATHETER, STRAIGHT	10/BOX
DSP136C	13F X 15CM DUO SPLIT® CATHETER, STRAIGHT	10/BOX
DSP138C	13F X 20CM DUO SPLIT® CATHETER, STRAIGHT	10/BOX
DSP139C	13F X 24CM DUO SPLIT® CATHETER, STRAIGHT	10/BOX
DSP130C	13F X 30CM DUO SPLIT® CATHETER, STRAIGHT	10/BOX
DSP134PC	13F X 12CM DUO SPLIT® CATHETER, PRE-CURVED	10/BOX
DSP136PC	13F X 15CM DUO SPLIT® CATHETER, PRE-CURVED	10/BOX
DSP138PC	13F X 20CM DUO SPLIT® CATHETER, PRE-CURVED	10/BOX
DSP139PC	13F X 24CM DUO SPLIT® CATHETER, PRE-CURVED	10/BOX
DDSP136IJC	13F X 15CM DUO SPLIT® CATHETER, RAULERSON IJ	10/BOX
DDSP138IJC	13F X 20CM DUO SPLIT® CATHETER, RAULERSON IJ	10/BOX
DDSP139IJC	13F X 24CM DUO SPLIT® CATHETER, RAULERSON IJ	10/BOX

CATHETER ONLY CONTENTS: (1) Catheter with Dual Stylets (2) Dilators (2) End Caps

BASIC SET

DDSP134S	13F X 12CM DUO SPLIT® CATHETER, STRAIGHT	5/BOX
DDSP136S	13F X 15CM DUO SPLIT® CATHETER, STRAIGHT	5/BOX
DDSP138S	13F X 20CM DUO SPLIT® CATHETER, STRAIGHT	5/BOX
DDSP139S	13F X 24CM DUO SPLIT® CATHETER, STRAIGHT	5/BOX
DDSP130S	13F X 30CM DUO SPLIT® CATHETER, STRAIGHT	5/BOX
DDSP134PS	13F X 12CM DUO SPLIT® CATHETER, PRE-CURVED	5/BOX
DDSP136PS	13F X 15CM DUO SPLIT® CATHETER, PRE-CURVED	5/BOX
DDSP138PS	13F X 20CM DUO SPLIT® CATHETER, PRE-CURVED	5/BOX
DDSP139PS	13F X 24CM DUO SPLIT® CATHETER, PRE-CURVED	5/BOX
DDSP134IJS	13F X 12CM DUO SPLIT® CATHETER, RAULERSON IJ	5/BOX
DDSP136IJS	13F X 15CM DUO SPLIT® CATHETER, RAULERSON IJ	5/BOX
DDSP138IJS	13F X 20CM DUO SPLIT® CATHETER, RAULERSON IJ	5/BOX
DDSP139IJS	13F X 24CM DUO SPLIT® CATHETER, RAULERSON IJ	5/BOX

BASIC SET CONTENTS: (1) Catheter with Dual Stylets (1) Introducer Needle (1) Guidewire (2) Dilators (1) Adhesive Wound Dressing (2) End Caps

FULL TRAY

DSP136T	13F X 15CM DUO SPLIT® CATHETER, STRAIGHT	5/BOX
DSP138T	13F X 20CM DUO SPLIT® CATHETER, STRAIGHT	5/BOX
DSP139T	13F X 24CM DUO SPLIT® CATHETER, STRAIGHT	5/BOX
DSP130T	13F X 30CM DUO SPLIT® CATHETER, STRAIGHT	5/BOX
DSP134PT	13F X 12CM DUO SPLIT® CATHETER, PRE-CURVED	5/BOX
DSP136PT	13F X 15CM DUO SPLIT® CATHETER, PRE-CURVED	5/BOX
DSP138PT	13F X 20CM DUO SPLIT® CATHETER, PRE-CURVED	5/BOX
DSP139PT	13F X 24CM DUO SPLIT® CATHETER, PRE-CURVED	5/BOX
DSP134IJT	13F X 12CM DUO SPLIT® CATHETER, RAULERSON IJ	5/BOX
DSP136IJT	13F X 15CM DUO SPLIT® CATHETER, RAULERSON IJ	5/BOX
DSP138IJT	13F X 20CM DUO SPLIT® CATHETER, RAULERSON IJ	5/BOX
DSP139IJT	13F X 24CM DUO SPLIT® CATHETER, RAULERSON IJ	5/BOX

FULL TRAY CONTENTS: (1) Catheter with Dual Stylets (1) Introducer Needle (1) Guidewire (2) Dilators (1) Adhesive Wound Dressing (2) End Caps (2) Chloraprep (1) 5ml Lidocaine (2) Pre-Filled Saline Syringes (1) 5cc Syringe (3) 4x4 Gauze (1) Surgical Gloves (1) Scalpel (3) Safety Hypodermic Needles (1) Filter Straw (1) 2-0 Silk Sutures (1) Fenestrated Drape (1) Hemostat (1) Lidocaine/Saline Sticker Sheet (1) Needle Stick Pad

IMPORTANT RISK INFORMATION

Indications for Use: The Medcomp® Duo-Split® Double Lumen Catheter is designed for acute hemodialysis and apheresis. It may be inserted percutaneously and is ideally placed in the internal jugular vein. Although this catheter may be inserted into the subclavian vein or femoral vein, the internal jugular vein is the preferred site.

Contraindications: This catheter is not intended for any use other than that which is indicated. Do not insert catheter in thrombosed vessels.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.