

DUO-SPLIT[®]

Short-Term Hemodialysis Catheter



FEATURES & BENEFITS

LUMEN DESIGN

Double "D" Internal



TIP DESIGN

Split Tip



AVAILABLE IN

13F



CONFIGURATIONS

Straight



Pre-Curved



Raulerson IJ



DUO-SPLIT®

Short-Term Hemodialysis Catheter

BASIC SET

| | | |
|-----------|--|-------|
| DDSP134S | 13F X 12CM DUO SPLIT® CATHETER, STRAIGHT | 5/BOX |
| DDSP136S | 13F X 15CM DUO SPLIT® CATHETER, STRAIGHT | 5/BOX |
| DDSP138S | 13F X 20CM DUO SPLIT® CATHETER, STRAIGHT | 5/BOX |
| DDSP139S | 13F X 24CM DUO SPLIT® CATHETER, STRAIGHT | 5/BOX |
| DDSP130S | 13F X 30CM DUO SPLIT® CATHETER, STRAIGHT | 5/BOX |
| DDSP134PS | 13F X 12CM DUO SPLIT® CATHETER, PRE-CURVED | 5/BOX |
| DDSP136PS | 13F X 15CM DUO SPLIT® CATHETER, PRE-CURVED | 5/BOX |
| DDSP138PS | 13F X 20CM DUO SPLIT® CATHETER, PRE-CURVED | 5/BOX |
| DDSP139PS | 13F X 24CM DUO SPLIT® CATHETER, PRE-CURVED | 5/BOX |
| DDSP134JS | 13F X 12CM DUO SPLIT® CATHETER, RAULERSON IJ | 5/BOX |
| DDSP136JS | 13F X 15CM DUO SPLIT® CATHETER, RAULERSON IJ | 5/BOX |
| DDSP138JS | 13F X 20CM DUO SPLIT® CATHETER, RAULERSON IJ | 5/BOX |
| DDSP139JS | 13F X 24CM DUO SPLIT® CATHETER, RAULERSON IJ | 5/BOX |

BASIC SET CONTENTS: (1) Catheter with Dual Stylets (1) Introducer Needle (1) Guidewire (2) Dilators (1) Adhesive Wound Dressing (2) End Caps

FULL TRAY

| | | |
|----------|--|-------|
| DSP136T | 13F X 15CM DUO SPLIT® CATHETER, STRAIGHT | 5/BOX |
| DSP138T | 13F X 20CM DUO SPLIT® CATHETER, STRAIGHT | 5/BOX |
| DSP139T | 13F X 24CM DUO SPLIT® CATHETER, STRAIGHT | 5/BOX |
| DSP130T | 13F X 30CM DUO SPLIT® CATHETER, STRAIGHT | 5/BOX |
| DSP134PT | 13F X 12CM DUO SPLIT® CATHETER, PRE-CURVED | 5/BOX |
| DSP136PT | 13F X 15CM DUO SPLIT® CATHETER, PRE-CURVED | 5/BOX |
| DSP138PT | 13F X 20CM DUO SPLIT® CATHETER, PRE-CURVED | 5/BOX |
| DSP139PT | 13F X 24CM DUO SPLIT® CATHETER, PRE-CURVED | 5/BOX |
| DSP134JT | 13F X 12CM DUO SPLIT® CATHETER, RAULERSON IJ | 5/BOX |
| DSP136JT | 13F X 15CM DUO SPLIT® CATHETER, RAULERSON IJ | 5/BOX |
| DSP138JT | 13F X 20CM DUO SPLIT® CATHETER, RAULERSON IJ | 5/BOX |
| DSP139JT | 13F X 24CM DUO SPLIT® CATHETER, RAULERSON IJ | 5/BOX |

FULL TRAY CONTENTS: (1) Catheter with Dual Stylets (1) Introducer Needle (1) Guidewire (2) Dilators (1) Adhesive Wound Dressing (2) End Caps (2) Chloraprep (1) 5ml Lidocaine (2) Pre-Filled Saline Syringes (1) 5cc Syringe (3) 4x4 Gauze (1) Surgical Gloves (1) Scalpel (3) Safety Hypodermic Needles (1) Filter Straw (1) 2-0 Silk Sutures (1) Fenestrated Drape (1) Hemostat (1) Lidocaine/Saline Sticker Sheet (1) Needle Stick Pad

IMPORTANT RISK INFORMATION

Indications for Use: The Medcomp® Duo-Split® Double Lumen Catheter is designed for acute hemodialysis and apheresis. It may be inserted percutaneously and is ideally placed in the internal jugular vein. Although this catheter may be inserted into the subclavian vein or femoral vein, the internal jugular vein is the preferred site.

Contraindications: This catheter is not intended for any use other than that which is indicated. Do not insert catheter in thrombosed vessels.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.