



HEMO-FLOW®/ HEMO-FLOW® XF

Indications for Use: The Hemo-Flow® and Hemo-Flow® XF Catheters are indicated for use in attaining Long-Term vascular access for Hemodialysis and Apheresis. It may be inserted percutaneously and is primarily placed in the internal jugular vein of an adult patient. Alternate insertion sites include subclavian vein as required. Catheters greater than 40cm are intended for femoral vein insertion.

Contraindications: These catheters are intended for Long-Term vascular access only and should not be used for any purpose other than indicated in these instructions.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.



Polyurethane
Material

Step Tip
Tip Design

14.5F
French Size(s)

**Straight,
Pre-Curved**
Configurations

Antegrade
Insertion Type



DOUBLE "D" INTERNAL LUMEN DESIGN



Radiopaque Polyurethane



ANGLED ARTERIAL TIP



ARTERIAL SIDEHOLES



BASIC SET		
HFS24E.	14.5F X 24CM HEMO-FLOW®, STRAIGHT	5/BOX
HFS28E.	14.5F X 28CM HEMO-FLOW®, STRAIGHT	5/BOX
HFS32E.	14.5F X 32CM HEMO-FLOW®, STRAIGHT	5/BOX
HFS36E.	14.5F X 36CM HEMO-FLOW®, STRAIGHT	5/BOX
HFS40E.	14.5F X 40CM HEMO-FLOW®, STRAIGHT	5/BOX
HFS55E.	14.5F X 55CM HEMO-FLOW®, STRAIGHT	5/BOX
HFS24PCE.	14.5F X 24CM PC HEMO-FLOW®, PRE-CURVED	5/BOX
HFS28PCE.	14.5F X 28CM PC HEMO-FLOW®, PRE-CURVED	5/BOX
HFS32PCE.	14.5F X 32CM PC HEMO-FLOW®, PRE-CURVED	5/BOX
HFS36PCE.	14.5F X 36CM PC HEMO-FLOW®, PRE-CURVED	5/BOX

- BASIC SET CONTAINS:**
- (1) Catheter
 - (1) Introducer Needle
 - (1) Guidewire
 - (2) Dilators
 - (1) Tuner
 - (1) Valved Peelable Introducer
 - (2) End Caps
 - (1) Scalpel

Refer to the Table of Contents for Important Risk Information regarding this device.