



LONG TERM PEDIATRIC SPLIT CATH® & SPLIT CATH® III

Polyurethane
Material

Split Tip
Tip Design


10F, 14F, 16F
French Sizes

**Straight,
Pre-Curved**
Configurations

Antegrade
Insertion Type



DOUBLE "D" INTERNAL
LUMEN DESIGN



ARTERIOVENOUS STYLET WEAVE



HOLLOW, FAST TRACK STYLET



SPLITTABLE LUMENS



Refer to the Table of Contents for Important Risk Information regarding this device.