




<b>Silicone</b> Material	<b>Straight, Coiled</b> Tip Designs	<b>15F</b> French Size	<b>Double</b> Cuff
-----------------------------	--	---------------------------	-----------------------




15F CATHETER ONLY STRAIGHT		
MC20VS43LC	15F X 43CM V-SERIES PERITONEAL DIALYSIS CATHETER, STRAIGHT W/ DOUBLE CUFF (LEFT)	5/BOX
MC20VS43RC	15F X 43CM V-SERIES PERITONEAL DIALYSIS CATHETER, STRAIGHT W/ DOUBLE CUFF (RIGHT)	5/BOX

15F CATHETER ONLY COILED		
MC20VC57LC	15F X 57CM V-SERIES PERITONEAL DIALYSIS CATHETER, COILED W/ DOUBLE CUFF (LEFT)	5/BOX
MC20VC63LC	15F X 63CM V-SERIES PERITONEAL DIALYSIS CATHETER, COILED W/ DOUBLE CUFF (LEFT)	5/BOX
MC20VC63RC	15F X 63CM V-SERIES PERITONEAL DIALYSIS CATHETER, COILED W/ DOUBLE CUFF (RIGHT)	5/BOX

15F BASIC SET STRAIGHT		
MC20VS43LS	15F X 43CM V-SERIES PERITONEAL DIALYSIS CATHETER, STRAIGHT W/ DOUBLE CUFF (LEFT)	5/BOX
MC20VS43RS	15F X 43CM V-SERIES PERITONEAL DIALYSIS CATHETER, STRAIGHT W/ DOUBLE CUFF (RIGHT)	5/BOX

15F BASIC SET COILED		
MC20VC57LS	15F X 57CM V-SERIES PERITONEAL DIALYSIS CATHETER, COILED W/ DOUBLE CUFF (LEFT)	5/BOX
MC20VC63LS	15F X 63CM V-SERIES PERITONEAL DIALYSIS CATHETER, COILED W/ DOUBLE CUFF (LEFT)	5/BOX
MC20VC63RS	15F X 63CM V-SERIES PERITONEAL DIALYSIS CATHETER, COILED W/ DOUBLE CUFF (RIGHT)	5/BOX

- |                                |                            |
|--------------------------------|----------------------------|
| <b>CATHETER ONLY CONTENTS:</b> | <b>BASIC SET CONTENTS:</b> |
| (1) Catheter                   | (1) Catheter               |
| (1) Female Luer                | (1) Female Luer            |
| (1) Clamp                      | (1) Clamp                  |
| (1) End Cap                    | (1) Introducer Needle      |
|                                | (1) Guidewire              |
|                                | (1) Peelable Introducer    |
|                                | (1) Tunneler               |
|                                | (3) 4x4 Gauze              |
|                                | (1) 10cc Syringe           |
|                                | (1) End Cap                |



Refer to the Table of Contents for Important Risk Information regarding this device.