

X-SERIES


Indications for Use: The Peritoneal Catheter is indicated for chronic peritoneal dialysis. The extended length of the “X” Series Catheter makes it especially applicable for peritoneal dialysis patients when it is necessary to locate the skin exit site remote from the usual lower abdominal region. The catheter may be particularly indicated in patients with obesity, floppy abdominal skin folds, urinary or fecal incontinence, chronic yeast intertrigo, intestinal stomas, or in patients who desire to take deep tub baths.

Contraindications: This catheter is contraindicated in patients who are not suitable candidates for peritoneal dialysis therapy. An upper chest (presteral) exit-site location of the X-Series Catheter may be contraindicated in patients who have had breast augmentation or breast reconstruction, or have a tracheostomy. However, this group of patients may be suitable candidates for an upper abdominal exit-site location of the X-Series Catheter.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.



Silicone Material	Coiled Tip Design	15F French Size
-----------------------------	-----------------------------	---------------------------



COILED TIP

LARGE SIDE HOLES

TITANIUM BARBED CONNECTOR



CATHETER ONLY		
MC20XC	15F X 64CM X-SERIES, COILED	5/BOX

BASIC SET		
MC20XS	15F X 64CM X-SERIES, COILED	5/BOX

ACCESSORIES		
MP217	17F VASCU-SHEATH® INTRODUCER SET	5/BOX
MC50	MALE LUER LOCK END CAP	100/BOX
MBCA	BARBED LUER LOCK CATHETER ADAPTOR	20/BOX
MCS	CATHETER STYLET	5/BOX
MTS	TUNNELING STYLET	10/BOX
MTC	TITANIUM LUER ADAPTOR	5/BOX
TBC	TITANIUM BARBED CONNECTOR	50/BOX
TRX	X-SERIES TUNNELING ROD	5/BOX
FSX	X-SERIES FALLER STYLET	5/BOX
40501	PREOPERATIVE MARKING STENCIL PACK	1/BOX

CATHETER ONLY CONTAINS:

- (1) Catheter
- (1) Titanium Barbed Connector
- (1) Female Luer
- (1) Clamp
- (1) End Cap

BASIC SET CONTAINS:

- (1) Catheter
- (1) Titanium Barbed Connector
- (1) Female Luer
- (1) Clamp
- (1) End Cap
- (1) Introducer Needle
- (1) Guidewire
- (1) Scalpel
- (1) Tunneler
- (1) 10cc Syringe
- (6) 4" x 4" Gauze
- (1) Tape Measure

Refer to the Table of Contents for Important Risk Information regarding this device.