



PRO-PICC®

Indications for Use: The Pro-PICC® CT catheter is indicated for short or long term peripheral access to the central venous system for intravenous therapy and power injection of contrast media, and allows for central venous pressure monitoring. For blood sampling, infusion, or therapies, use a 4F or larger catheter. The maximum recommended infusion rate varies by catheter French size and is printed on the catheter.

Contraindications: The presence of device related infection, bacteremia, or septicemia is known or suspected. The patient's body size is insufficient to accommodate the size of the implanted device. The patient is known or is suspected to be allergic to materials contained in the device. There has been past irradiation of prospective insertion site. There have been previous episodes of venous thrombosis or vascular surgical procedures at the prospective placement site. There are local tissue factors that may prevent proper device stabilization and/or access.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.



Single,
Double,
Triple
Lumens

55cm, 60cm
Lengths


3F, 4F, 5F, 6F
French Sizes

5cc/sec
@300psi*
Injection Rating

PICCs



POWER INJECTABLE
Capable of CT injections
up to 5cc/sec at 300 psi for
diagnostic imaging.*



*Except 3F Pro-PICC®



Refer to the Table of Contents for Important Risk Information regarding this device.



CATHETER ONLY		
MR17033100	3F X 55CM PRO-PICC®, SINGLE LUMEN	EACH
MR17034100	4F X 55CM PRO-PICC®, SINGLE LUMEN	EACH
MR17035100	5F X 60CM PRO-PICC®, SINGLE LUMEN	EACH
MR17035200	5F X 55CM PRO-PICC®, DOUBLE LUMEN	EACH

BASIC SET		
MR17033101	3F X 55CM PRO-PICC®, SINGLE LUMEN	EACH
MR17034101	4F X 55CM PRO-PICC®, SINGLE LUMEN	EACH
MR17035101	5F X 60CM PRO-PICC®, SINGLE LUMEN	EACH
MR17035201	5F X 55CM PRO-PICC®, DOUBLE LUMEN	EACH
MR17036201	6F X 60CM PRO-PICC®, DOUBLE LUMEN	EACH
MR17036301	6F X 60CM PRO-PICC®, TRIPLE LUMEN	EACH

BASIC SET WITH FLUSHABLE STYLET		
MR17034111	4F X 55CM PRO-PICC®, SINGLE LUMEN	EACH
MR17035111	5F X 60CM PRO-PICC®, SINGLE LUMEN	EACH
MR17035211	5F X 55CM PRO-PICC®, DOUBLE LUMEN	EACH
MR17036311	6F X 60CM PRO-PICC®, TRIPLE LUMEN	EACH

LONG WIRE SET		
MR17033102	3F X 55CM PRO-PICC®, SINGLE LUMEN	EACH
MR17034102	4F X 55CM PRO-PICC®, SINGLE LUMEN	EACH
MR17035102	5F X 60CM PRO-PICC®, SINGLE LUMEN	EACH
MR17035202	5F X 55CM PRO-PICC®, DOUBLE LUMEN	EACH
MR17036202	6F X 60CM PRO-PICC®, DOUBLE LUMEN	EACH
MR17036302	6F X 60CM PRO-PICC®, TRIPLE LUMEN	EACH

NURSING KIT		
MR17033105	3F X 55CM PRO-PICC®, SINGLE LUMEN	EACH
MR17034105	4F X 55CM PRO-PICC®, SINGLE LUMEN	EACH
MR17035105	5F X 60CM PRO-PICC®, SINGLE LUMEN	EACH
MR17035205	5F X 55CM PRO-PICC®, DOUBLE LUMEN	EACH
MR17036205	6F X 60CM PRO-PICC®, DOUBLE LUMEN	EACH
MR17036305	6F X 60CM PRO-PICC®, TRIPLE LUMEN	EACH

CUT DOWN VERSION SAFETY		
MRCTP41024	4F X 55CM PRO-PICC®, SINGLE LUMEN	EACH
MRCTP51024	5F X 60CM PRO-PICC®, SINGLE LUMEN	EACH
MRCTP52024	5F X 55CM PRO-PICC®, DOUBLE LUMEN	EACH
MRCTP62024	6F X 60CM PRO-PICC®, DOUBLE LUMEN	EACH
MRCTP63024	6F X 60CM PRO-PICC®, TRIPLE LUMEN	EACH

CUT DOWN VERSION NON SAFETY		
MRCTP41028	4F X 55CM PRO-PICC®, SINGLE LUMEN	EACH
MRCTP51028	5F X 60CM PRO-PICC®, SINGLE LUMEN	EACH
MRCTP52028	5F X 55CM PRO-PICC®, DOUBLE LUMEN	EACH
MRCTP62028	6F X 60CM PRO-PICC®, DOUBLE LUMEN	EACH

CATHETER ONLY CONTAINS:
 (1) PICC w/ Stylet
 (1|2) Needle Free Valve(s)
 (1) Catheter Securement Device
 (1) Patient Information Pack
 (1) Patient Chart Sticker
 (1) Patient ID Card

BASIC SET CONTAINS:
 (1) PICC w/ Stylet
 (1) Peelable Sheath Introducer
 (1) Mini Scalpel
 (1) Guidewire
 (1) 10cc Syringe
 (1) Tape Measure
 (1) Introducer Needle
 (1|2|3) Needle Free Valve(s)
 (1) Catheter Securement Device
 (1) Patient Information Pack
 (1) Patient Chart Sticker
 (1) Patient ID Card

BASIC SET w/ FLUSHABLE STYLET CONTAINS:
 (1) PICC w/ Flushable Stylet
 (1) Peelable Sheath Introducer
 (1) Mini Scalpel
 (1) Guidewire
 (1) 10cc Syringe
 (1) Tape Measure
 (1) Introducer Needle
 (1|2|3) Needle Free Valve(s)
 (1) Catheter Securement Device
 (1) Patient Information Pack
 (1) Patient Chart Sticker
 (1) Patient ID Card

LONG WIRE SET CONTAINS:
 (1) PICC
 (1) Peelable Sheath Introducer
 (1) Mini Scalpel
 (1) Guidewire
 (1) 10cc Syringe
 (1) Tape Measure
 (1) Introducer Needle
 (1|2|3) Needle Free Valve(s)
 (1) Catheter Securement Device
 (1) Patient Information Pack
 (1) Patient Chart Sticker
 (1) Patient ID Card

NURSING KIT CONTAINS:
 (1) PICC w/ Stylet
 (1) Peelable Sheath Introducer
 (1) Safety Scalpel
 (1) Guidewire
 (1) 10cc Syringe
 (1) Tape Measure
 (1) Safety Introducer Needle
 (1|2|3) Needle Free Valve(s)
 (1) Latex-Free Tourniquet
 (1) Catheter Securement Device
 (1) Patient Information Pack
 (1) Patient Chart Sticker
 (1) Patient ID Card

CUT DOWN VERSION SAFETY CONTAINS:
 (1) PICC w/ Stylet
 (1) Peelable Sheath Introducer
 (1) Safety Introducer Needle
 (1) Guidewire
 (1) Safety Scalpel
 (1) Tape Measure
 (1) Catheter Securement Device

CUT DOWN VERSION NON SAFETY CONTAINS:
 (1) PICC w/ Stylet
 (1) Peelable Sheath Introducer
 (1) Introducer Needle
 (1) Guidewire
 (1) Mini Scalpel
 (1) Tape Measure
 (1) Catheter Securement Device

