



VASCU-PICC® / VASCU-PICC® TL

**Indications for Use:** The Peripherally Inserted Central Vein Access Catheters are designed for Long or Short- Term peripheral access to the central venous system for intravenous therapy and blood sampling, and allows for central venous pressure monitoring. This catheter may be inserted via the basilic, cephalic, or median cubital vein.

**Contraindications:** This catheter is not intended for any use other than that which is indicated. Do not implant catheter in thrombosed vessels. The presence of skin related problems around the insertion site (infection, phlebitis, scars, etc.). The presence of device related bacteremia or septicemia. History of mastectomy on insertion side. Previous history of venous/subclavian thrombosis or vascular surgical procedures at insertion site. Fever of unknown origin. The patient's body size is insufficient to accommodate the size of the implanted device. The patient is known or is suspected to be allergic to materials contained in the device. Past irradiation of prospective insertion site. Local tissue factors will prevent proper device stabilization and/or access.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.



<b>Single, Double</b> Lumens	<b>60cm</b> Length	<b>3F, 4F, 5F</b> French Sizes
-------------------------------------	-----------------------	-----------------------------------



**TAPERLESS LUMEN DESIGN**



**DEPTH MARKS**



CATHETER ONLY		
MR81013100	3F X 60CM VASCU-PICC® TL, SINGLE LUMEN	EACH
MR81014100	4F X 60CM VASCU-PICC® TL, SINGLE LUMEN	EACH
MR81014200	4F X 60CM VASCU-PICC® TL, DOUBLE LUMEN	EACH

BASIC SET		
MR81013101	3F X 60CM VASCU-PICC® TL, SINGLE LUMEN	EACH
MR81014101	4F X 60CM VASCU-PICC® TL, SINGLE LUMEN	EACH
MR81015101	5F X 60CM VASCU-PICC® TL, SINGLE LUMEN	EACH
MR81014201	4F X 60CM VASCU-PICC® TL, DOUBLE LUMEN	EACH
MR81015201	5F X 60CM VASCU-PICC® TL, DOUBLE LUMEN	EACH

BASIC NURSING		
MR81013105	3F X 60CM VASCU-PICC® TL, SINGLE LUMEN	EACH
MR81014105	4F X 60CM VASCU-PICC® TL, SINGLE LUMEN	EACH
MR81015105	5F X 60CM VASCU-PICC® TL, SINGLE LUMEN	EACH
MR81014205	4F X 60CM VASCU-PICC® TL, DOUBLE LUMEN	EACH
MR81015205	5F X 60CM VASCU-PICC® TL, DOUBLE LUMEN	EACH

**CATHETER ONLY CONTAINS:**  
 (1) Vascu-PICC® TL w/ Stylet  
 (1|2) Needle Free Valve(s)  
 (1) Catheter Securement Device  
 (1) Patient Information Pack  
 (1) Patient Chart Sticker

**BASIC SET CONTAINS:**  
 (1) Vascu-PICC® TL w/ Stylet  
 (1) Introducer Needle  
 (1) Peelable Sheath Introducer  
 (1) 10cc Syringe  
 (1) Guidewire  
 (1) Mini Scalpel  
 (1|2) Needle Free Valve(s)  
 (1) Tape Measure  
 (1) Catheter Securement Device  
 (1) Patient Information Pack  
 (1) Patient Chart Sticker

**BASIC NURSING CONTAINS:**  
 (1) Vascu-PICC® TL w/ Stylet  
 (1) Safety Introducer Needle  
 (1) Tape Measure  
 (1) Peelable Sheath Introducer  
 (1) 10cc Syringe  
 (1) Guidewire  
 (1) Safety Scalpel  
 (1|2) Needle Free Valve(s)  
 (1) Latex Free Tourniquet  
 (1) Catheter Securement Device  
 (1) Patient Information Pack  
 (1) Patient Chart Sticker

Refer to the Table of Contents for Important Risk Information regarding this device.