



VASCU-PICC® / VASCU-PICC® TL

Indications for Use: The Peripherally Inserted Central Vein Access Catheters are designed for Long or Short- Term peripheral access to the central venous system for intravenous therapy and blood sampling, and allows for central venous pressure monitoring. This catheter may be inserted via the basilic, cephalic, or median cubital vein.

Contraindications: This catheter is not intended for any use other than that which is indicated. Do not implant catheter in thrombosed vessels. The presence of skin related problems around the insertion site (infection, phlebitis, scars, etc.). The presence of device related bacteremia or septicemia. History of mastectomy on insertion side. Previous history of venous/subclavian thrombosis or vascular surgical procedures at insertion site. Fever of unknown origin. The patient's body size is insufficient to accommodate the size of the implanted device. The patient is known or is suspected to be allergic to materials contained in the device. Past irradiation of prospective insertion site. Local tissue factors will prevent proper device stabilization and/or access.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.

Single, Double Lumens	60cm Length	3F, 4F, 5F, 6F French Sizes
-------------------------------------	-----------------------	---------------------------------------


REVERSE TAPER LUMEN DESIGN

DEPTH MARKS


CATHETER ONLY		
MR17013100	3F X 60CM VASCU-PICC®, SINGLE LUMEN	EACH
MR17014100	4F X 60CM VASCU-PICC®, SINGLE LUMEN	EACH
MR17014200	4F X 60CM VASCU-PICC®, DOUBLE LUMEN	EACH
MR17015200	5F X 60CM VASCU-PICC®, DOUBLE LUMEN	EACH
MR17016200	6F X 60CM VASCU-PICC®, DOUBLE LUMEN	EACH

BASIC SET		
MR17013101	3F X 60CM VASCU-PICC®, SINGLE LUMEN	EACH
MR17014101	4F X 60CM VASCU-PICC®, SINGLE LUMEN	EACH
MR17015101	5F X 60CM VASCU-PICC®, SINGLE LUMEN	EACH
MR17014201	4F X 60CM VASCU-PICC®, DOUBLE LUMEN	EACH
MR17015201	5F X 60CM VASCU-PICC®, DOUBLE LUMEN	EACH
MR17016201	6F X 60CM VASCU-PICC®, DOUBLE LUMEN	EACH

LONG WIRE SET		
MR17014102	4F X 60CM VASCU-PICC®, SINGLE LUMEN	EACH
MR17015202	5F X 60CM VASCU-PICC®, DOUBLE LUMEN	EACH

NURSING KIT		
MR17013105	3F X 60CM VASCU-PICC®, SINGLE LUMEN	EACH
MR17014105	4F X 60CM VASCU-PICC®, SINGLE LUMEN	EACH
MR17015105	5F X 60CM VASCU-PICC®, SINGLE LUMEN	EACH
MR17014205	4F X 60CM VASCU-PICC®, DOUBLE LUMEN	EACH
MR17015205	5F X 60CM VASCU-PICC®, DOUBLE LUMEN	EACH

OTN SET		
VP3S-MNS	3F X 60CM SINGLE VASCU-PICC® W/ OTN INTRODUCER SET	EACH
VP5S-MNS	5F X 60CM SINGLE VASCU-PICC® W/ OTN INTRODUCER SET	EACH
VP4D-MNS	4F X 60CM DUAL VASCU-PICC® W/ OTN INTRODUCER SET	EACH
VP5D-MNS	5F X 60CM DUAL VASCU-PICC® W/ OTN INTRODUCER SET	EACH

CATHETER ONLY CONTAINS:
 (1) PICC w/ Stylet
 (1|2) Needle Free Valve(s)
 (1) Catheter Securement Device
 (1) Patient Information Pack
 (1) Patient Chart Sticker
 (1) Vascu-PICC® IFU

BASIC SET CONTAINS:
 (1) PICC w/ Stylet
 (1) Peelable Sheath Introducer
 (1) Mini Scalpel
 (1) Guidewire
 (1) 10cc Syringe
 (1) Tape Measure
 (1) Introducer Needle
 (1|2) Needle Free Valve(s)
 (1) Catheter Securement Device
 (1) Patient Information Pack
 (1) Patient Chart Sticker

LONG WIRE SET CONTAINS:
 (1) PICC
 (1) Peelable Sheath Introducer
 (1) Mini Scalpel
 (1) Guidewire
 (1) 10cc Syringe
 (1) Tape Measure
 (1) Introducer Needle
 (1|2) Needle Free Valve(s)
 (1) Catheter Securement Device
 (1) Patient Information Pack
 (1) Patient Chart Sticker

NURSING KIT CONTAINS:
 (1) PICC w/ Stylet
 (1) Peelable Sheath Introducer
 (1) Safety Scalpel
 (1) Guidewire
 (1) 10cc Syringe
 (1) Tape Measure
 (1) Safety Introducer Needle
 (1|2) Needle Free Valve(s)
 (1) Latex-Free Tourniquet
 (1) Catheter Securement Device
 (1) Patient Information Pack
 (1) Patient Chart Sticker

OTN SET CONTAINS:
 (1) Vascu-PICC® TLS
 (1) Peelable Sheath Introducer
 (1) Guidewire
 (1) Tape Measure
 (1) 10cc Syringe
 (1|2) Needle Free Valve(s)
 (1) Catheter Securement Device
 (1) Patient Chart Sticker
 (1) Patient Information Pack

Refer to the Table of Contents for Important Risk Information regarding this device.