



DIGNITY® CT/ DIGNITY® MINI/ DIGNITY® DMP

Indications for Use: The CT Power Injectable Implantable Infusion Ports are indicated for patient therapies requiring repeated access to the vascular system. The port system can be used for infusion of medications, I.V. fluids, parenteral nutrition solutions, blood products, and for the withdrawal of blood samples. When used with a power injectable needle, the Power Injectable Implantable Infusion Port device is indicated for power injection of contrast media. For power injection of contrast media, the maximum recommended infusion rate is 5ml/s with a 19 or 20 gauge non-coring power injectable needle. The maximum recommended infusion rate is 2ml/s with a 22 gauge non-coring power injectable needle.

Contraindications: This device is contraindicated for catheter insertion in the subclavian vein medial to the border of the first rib, an area which is associated with higher rates of pinch-off. The device is also contraindicated: When the presence of device related infection, bacteremia, or septicemia is known or suspected. When the patient's body size is insufficient for the size of the implanted device. When the patient is known or is suspected to be allergic to materials contained in the device. If severe chronic obstructive lung disease exists. If the prospective insertion site has been previously irradiated. If the prospective placement site has previously suffered episodes of venous thrombosis or vascular surgical procedures. If local tissue factors will prevent proper device stabilization and/or access.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.



Plastic
Port Material

**Mini,
Low Profile,
Mid Sized**
Configurations

**Polyurethane,
Silicone**
Catheter Materials

**5cc/sec
@ 325psi**
Injection Rating

**Attachable,
Pre-Attached**
Catheter
Configurations

**5F, 6.6F,
8F, 9.6F**
French Sizes

Ports



LOCKING MECHANISM



PATIENT-CENTRIC DESIGN

SCULPTED SIDES



CT UNDER X-RAY



Refer to the Table of Contents for Important Risk Information regarding this device.

DIGNITY® LOW PROFILE			SILICONE FILLED SUTURE HOLES		
	WITHOUT MICRO INTRODUCER	WITH MICRO INTRODUCER	WITHOUT MICRO INTRODUCER	WITH MICRO INTRODUCER	
5F DIGNITY® LP SHAPED ATTACHABLE	MRCTI50041	MICTI50041	MRCTI5004S	MICTI5004S	5/BOX
5F DIGNITY® LP SHAPED PRE-ATTACHED	MRCTI50841	MICTI50841	MRCTI5084S	MICTI5084S	5/BOX
6.6F DIGNITY® LP SHAPED ATTACHABLE	MRCTI66041	MICTI66041	MRCTI6604S	MICTI6604S	5/BOX
6.6F DIGNITY® LP SHAPED PRE-ATTACHED	MRCTI66841	MICTI66841	MRCTI6684S	MICTI6684S	5/BOX
8F DIGNITY® LP SHAPED ATTACHABLE	MRCTI80041	MICTI80041	MRCTI8004S	MICTI8004S	5/BOX
8F DIGNITY® LP SHAPED PRE-ATTACHED	MRCTI80841	MICTI80841	MRCTI8084S	MICTI8084S	5/BOX
DIGNITY® MID SIZE					
5F DIGNITY® MS SHAPED ATTACHABLE	MRCTI50001	MICTI50001	MRCTI5000S	MICTI5000S	5/BOX
5F DIGNITY® MS SHAPED PRE-ATTACHED	MRCTI50801	MICTI50801	MRCTI5080S	MICTI5080S	5/BOX
6.6F DIGNITY® MS SHAPED ATTACHABLE	MRCTI66001	MICTI66001	MRCTI6600S	MICTI6600S	5/BOX
6.6F DIGNITY® MS SHAPED PRE-ATTACHED	MRCTI66801	MICTI66801	MRCTI6680S	MICTI6680S	5/BOX
8F DIGNITY® MS SHAPED ATTACHABLE	MRCTI80001	MICTI80001	MRCTI8000S	MICTI8000S	5/BOX
8F DIGNITY® MS SHAPED PRE-ATTACHED	MRCTI80801	MICTI80801	MRCTI8080S	MICTI8080S	5/BOX
9.6F DIGNITY® MS SHAPED ATTACHABLE	MRCTI96001	MICTI96001	MRCTI9600S	MICTI9600S	5/BOX
9.6F DIGNITY® MS SHAPED PRE-ATTACHED	MRCTI96801	MICTI96801	MRCTI9680S	MICTI9680S	5/BOX

 Silicone Catheters

SETS CONTAIN:

- (1) Catheter
- (1) CT Implantable Port
- (2) Catheter Locks
- (1) Scalpel
- (1) 18ga Introducer Needle
- (1) Vein Pick
- (1) .035" x 70cm "J" Marked Guidewire
- with Advancer
- (1) 10cc Syringes
- (1) Peelable Introducer
- (1) Tunneler
- (2) 22ga Huber Needle - (1 Straight, 1 Right Angle)
- (1) Blunt Tip Needle
- (1) Patient Chart Sticker
- (1) Patient Information Pack

SETS W/ MICRO-STICK®

- CONTAIN:
- (1) Catheter
 - (1) CT Implantable Port
 - (2) Catheter Locks
 - (1) Scalpel
 - (1) Vein Pick
 - (1) .035" x 70cm "J" Marked Guidewire
 - with Advancer
 - (1) 10cc Syringes
 - (1) Peelable Introducer
 - (1) Tunneler
 - (2) 22ga Huber Needle - (1 Straight, 1 Right Angle)
 - (1) Blunt Tip Needle
 - (1) Patient Chart Sticker
 - (1) Patient Information Pack

MICRO-STICK®

- (1) 5FR Coaxial Dilator
- (1) .018" x 45cm Guidewire with Radiopaque Tip
- (1) 21ga Introducer Needle with Echogenic Tip

