



DIGNITY® CT/ DIGNITY® MINI/ DIGNITY® DMP

Indications for Use: The CT Power Injectable Implantable Infusion Ports are indicated for patient therapies requiring repeated access to the vascular system. The port system can be used for infusion of medications, I.V. fluids, parenteral nutrition solutions, blood products, and for the withdrawal of blood samples. When used with a power injectable needle, the Power Injectable Implantable Infusion Port device is indicated for power injection of contrast media. For power injection of contrast media, the maximum recommended infusion rate is 5ml/s with a 19 or 20 gauge non-coring power injectable needle. The maximum recommended infusion rate is 2ml/s with a 22 gauge non-coring power injectable needle.

Contraindications: This device is contraindicated for catheter insertion in the subclavian vein medial to the border of the first rib, an area which is associated with higher rates of pinch-off. The device is also contraindicated: When the presence of device related infection, bacteremia, or septicemia is known or suspected. When the patient's body size is insufficient for the size of the implanted device. When the patient is known or is suspected to be allergic to materials contained in the device. If severe chronic obstructive lung disease exists. If the prospective insertion site has been previously irradiated. If the prospective placement site has previously suffered episodes of venous thrombosis or vascular surgical procedures. If local tissue factors will prevent proper device stabilization and/or access.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.


Plastic Port Material	Mini Configuration	Polyurethane Catheter Material	5cc/sec @ 325psi Injection Rating	Attachable, Pre-Attached Catheter Configurations	5F, 6.6F, 8F French Sizes
---------------------------------	------------------------------	--	---	--	-------------------------------------



REDUCED WIDTH AND SMALLER PROFILE

CHEST OR ARM PLACEMENT

PATIENT-CENTRIC DESIGN



DIGNITY® MINI	OPEN SUTURE HOLES		SILICONE FILLED SUTURE HOLES		
	WITHOUT MICRO INTRODUCER	WITH MICRO INTRODUCER	WITHOUT MICRO INTRODUCER	WITH MICRO INTRODUCER	
5F DIGNITY® MINI SHAPED ATTACHABLE	MRCTI50041M	MICTI50041M	MRCTI5004SM	MICTI5004SM	5/BOX
5F DIGNITY® MINI SHAPED PRE-ATTACHED	MRCTI50841M	MICTI50841M	MRCTI5084SM	MICTI5084SM	5/BOX
6.6F DIGNITY® MINI SHAPED ATTACHABLE	MRCTI66041M	MICTI66041M	MRCTI6604SM	MICTI6604SM	5/BOX
6.6F DIGNITY® MINI SHAPED PRE-ATTACHED	MRCTI66841M	MICTI66841M	MRCTI6684SM	MICTI6684SM	5/BOX
8F DIGNITY® MINI SHAPED ATTACHABLE	MRCTI80041M	MICTI80041M	MRCTI8004SM	MICTI8004SM	5/BOX
8F DIGNITY® MINI SHAPED PRE-ATTACHED	MRCTI80841M	MICTI80841M	MRCTI8084SM	MICTI8084SM	5/BOX

SET CONTAINS:

- (1) CT Implantable Port
- (1) Catheter
- (2) Catheter Locks
- (1) Scalpel
- (1) 18ga Introducer Needle
- (1) Vein Pick
- (1) 0.35" x 70cm "J" Marked Guidewire w/ Advancer
- (1) 10cc Syringe
- (1) Peelable Introducer
- (1) Tunneler
- (2) 22ga Huber Needles - (1 Straight, 1 Right Angle)
- (1) Blunt Tip Needle
- (1) Patient Chart Sticker
- (1) Patient Information Pack

SETS WITH MICRO-STICK® CONTAIN:

- (1) Catheter
- (1) CT Implantable Port
- (2) Catheter Locks
- (1) Scalpel
- (1) Vein Pick
- (1) .035" x 70cm "J" Marked Guidewire w/ Advancer
- (1) 10cc Syringe
- (1) Peelable Introducer
- (1) Tunneler
- (2) 22ga Huber Needles - (1 Straight, 1 Right Angle)
- (1) Blunt Tip Needle
- (1) Patient Chart Sticker
- (1) Patient Information Pack

MICRO-STICK®

- (1) 5FR Coaxial Dilator
- (1) .018" x 45cm Guidewire w/ Radiopaque Tip
- (1) 21ga Introducer Needle w/ Echogenic Tip

Refer to the Table of Contents for Important Risk Information regarding this device.