

DUO-FLOW® 400XL

Indications for Use: The Medcomp® Duo-Flow and Duo-Flow 400XL Double Lumen Catheters are designed for acute hemodialysis and apheresis. They may be inserted percutaneously and are ideally placed in the internal jugular vein. Although these catheters may be inserted into the subclavian or femoral vein, the internal jugular is the preferred site.

Contraindications: This catheter is intended for Short-Term vascular access only and should not be used for any purpose other than indicated in these instructions.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.



Polyurethane
Material

Tapered Tip
Tip Design

14F
French Size

**Straight,
Pre-Curved**
Configurations



COAXIAL INTERNAL LUMEN DESIGN



TAPERED TIP

360° ARTERIAL PORTS



BASIC SET		
DFXL146MTE	14F X 15CM DUO-FLOW® 400XL BASIC SET, STRAIGHT	10/BOX
DFXL148MTE	14F X 20CM DUO-FLOW® 400XL BASIC SET, STRAIGHT	10/BOX
DFXL149MTE	14F X 24CM DUO-FLOW® 400XL BASIC SET, STRAIGHT	10/BOX
DFXL144IJSE	14F X 12CM DUO-FLOW® 400XL BASIC SET, PRE-CURVED	10/BOX
DFXL146IJSE	14F X 15CM DUO-FLOW® 400XL BASIC SET, PRE-CURVED	10/BOX
DFXL148IJSE	14F X 20CM DUO-FLOW® 400XL BASIC SET, PRE-CURVED	10/BOX
DFXL149IJSE	14F X 24CM DUO-FLOW® 400XL BASIC SET, PRE-CURVED	10/BOX

**BASIC SET
CONTAINS:**
(1) Catheter
(1) Dilator
(1) Guidewire
(1) Introducer Needle
(2) End Caps

Refer to the Table of Contents for Important Risk Information regarding this device.